

Product datasheet for **RN215298**

Olr714 (NM_001000923) Rat Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Olr714 (NM_001000923) Rat Untagged Clone
Tag:	Tag Free
Symbol:	Olr714
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>RN215298 representing NM_001000923 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGGGCCAGAGCAAAATTGTCACAGAATTTGCTCTGCTGGGCTTCACTCAGGATCCTGCTGGGCAAAAAG
CTTTATTTGTTATGTTTTCACTGATGTACATTGCAACAATGGTGGGCAACCTGCTCATTGTAGGGACAGT
GATTGTCAGCCCTCCTTAGGCTCCCAATGTACTTCTTCCTTGCCTCTCTGTCACTAATGGATGCTGTT
TATTCTACTGCCATCTCACCCAAGTTGATTGTAGACTTGCTCCGTGAGAAAAAGACCATCTCCTTCAGAG
CTTGCCATAAGCCAGCTCTTTATAGAGCATTTATTTGGTGGGTTGACATTGTCATTCTGGTGGCAATGGC
CTATGATTGCTATGTGGCCATCTGTAAGCCACTACATTATTTGATCATCATGAATCGACGCGTATGTATC
CTCTTCTTGGTATGGCCTGGGAGGAGTTTTGCACATGCTTTGTTTCAAGTTCTTGCTGTATATAATC
TTCCTTCTGTGGTCCCAATATCATTGACCACTTTGGCTGTGACATATACCCGCTCTTGCTACTTGCATG
TACTGACACCTACCTGATTGGTCTCACTATCATTGCTAATGGTGGAGCCATGTGTATAGTCATCTTCATC
CTACTTCTACTTTCTTATGCAATAATCCTAAGATCTCTCCAGAATCACAGTCAGGAAGGGAGGCACAAAG
CTCTGTCCACCTGCAGCTCTCATATCACTGTGTTTTTCTCTTTTTTGTTCCTGTATTTAATGTATGT
TAGACCTGTTTCTAGTTTCTATTGACAAATCCATCACTGTAATTTATACTGTTGTACCCCCGTGTTG
AATCCTTTAATATATACTTTGAGAAATTCAGAGATGAAAAATTCTATGAAAAGCTCCTGCAGAAATAT
TAAGTCGAATAGACTAAGGCTTCTTTGCTGCTTATATTAT**TAA**

AG**CGGACCG**ACGCGTACGCGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCC
TGGATTACAAGGATGACGACGATAAGGTTTAA

Restriction Sites:	SgfI-RsrII
ACCN:	NM_001000923
Insert Size:	954 bp



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001000923.1</u> , <u>NP_001000923.1</u>
RefSeq Size:	954 bp
RefSeq ORF:	954 bp
Locus ID:	405239
Cytogenetics:	3q24
Gene Summary:	Olfactory receptors interact with odorant molecules in the nose, to initiate a neuronal response that triggers the perception of a smell. The olfactory receptor proteins are members of a large family of G-protein-coupled receptors (GPCR) arising from single coding-exon genes. Olfactory receptors share a 7-transmembrane domain structure with many neurotransmitter and hormone receptors and are responsible for the recognition and G protein-mediated transduction of odorant signals. The olfactory receptor gene family is the largest in the genome. The nomenclature assigned to the olfactory receptor genes and proteins for this organism is independent of other organisms. [provided by RefSeq, Jul 2008]