

Product datasheet for RN215337

Rsph1 (NM_001012176) Rat Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Rsph1 (NM_001012176) Rat Untagged Clone
Tag:	Tag Free
Symbol:	Rsph1
Synonyms:	Tsga2
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>RN215337 representing NM_001012176 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGTCGGACCTGGGCTCTGAGGAGTTAGAGGAGGAGGGAGAGAACGATCTTGGGAATACGAAGGAGAGC
GCAATGAACTGGGAGAGCGGCACGGACATGGGAAGGCACGGTTGCCAACGGGGACACATATGAAGGAAA
CTATGAGTCTGGAAAAAGACATGGCCAGGGCACCTACAAATTTAAGAATGGGGCCCGGTACGTTGGCGAC
TATGTTAAAAATAAGAAGCATGGACAAGGTACCTTTGTCTATCCAGATGGTTCCAGATATGAAGGGGAGT
GGGCTGACGACCAGAGGCATGGCCATGGCGTTTACTACTACGTCATAATGACACCTACACAGGGGAGTG
GTTCAATCATCAAAGGCACGGGCAAGGCACCTACTTCTACGGGAGACCGGCAGTAAGTACGTCGGTACC
TGGGTGCACGGGCAGCAGGAGGGTGCCGCCGAGCTCATCCACCTGAATCACAGGTACCAGGGCAAGTTCA
TAAACAAAACCCCTGTGGGCCCGGGAAAGTACGTATTTGACGTGGGATGCGAACAGCAGCGTGAATACCG
CCTCACAGATGTGGAAAGAGGTGAAGAGGAGGAGGAGGAAGACATTAGTGAGTATCGTCCAAAATGG
AAAGCTCTCAAAATCACAGAATTGGCCTTGTGGACTCCAACCCCTCCGAGGAGCATCCTCTGCTGAGG
GACAAGGGCAGGAAGAAGTCCAGGGAGCTACTGGTGTGGTGGCCCTCAGATGACATCCAGGCAGGAAGG
TTTTGAGGTCGACGCAGAGCCGAGGACGGGCGAAGAAGACATAGACTTCAGGCAGGAGAGCCAGGAGAAC
AGTTACGACTTAGACCAGGGAACCCGAACCTTTGATGAAGAACAGTTAGACCTCCAGGAGTAA

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites:	Sgfl-MluI
ACCN:	NM_001012176
Insert Size:	903 bp



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001012176.1, NP_001012176.1</u>
RefSeq Size:	1130 bp
RefSeq ORF:	903 bp
Locus ID:	361818
Cytogenetics:	20p12