

Product datasheet for RN215380

Serpinf1 (NM_177927) Rat Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Serpinf1 (NM_177927) Rat Untagged Clone
Tag:	Tag Free
Symbol:	Serpinf1
Synonyms:	Dmrs91; Pedf
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>RN215380 representing NM_177927 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGCAGACCCGTGGTACTCCTCTGGACTGGAGCCCTGCTTGGGCACGGCAGCAGCCAGAATGTCCTTG
ACAGCTCTCAGGATCCCCAGCCCTGACAGCACCGGGGAGCCCGTAGTGGAGGAGGATGACCCCTTCTT
CAAGGCCCCCGTGAACAAGTTGGCAGCAGCTGTTCCAACCTCGGCTACGATCTCTACCGCTGAGATCC
GGTGTCTCAACCGGCAACATCTGTGTCTCTCTCAGCGTGGCCACGGCCCTCTCGGCCCTTCC
TGGGAGCTGAACAGCGAACAGAGTCTGTCAATCACGGGCTCTACTACGACTTGATCAACAATCTGA
CATCCACAGCACCTACAAGGAGCTCCTTGCTCTGTTACTGCCCTGAGAAGAATTCAAGAGTGCCTCC
AGAATTGTGTTGAGAGGAACTTCGAGTAAAAATCCAGCTTTGTTGCTCCTCTGGAGAAGTCATATGGGA
CCAGGCCCCGAATCCTCACTGGCAACCCTCGCATAGACCTTCAGGAGATTAACAATGGGTGCAGGCCCA
GATGAAAGGGAAAATGCCCGGTCTACAAGGAAATGCCAGTGCCTCAGCATCCTCCTTGGCGTG
GCTTACTTCAAGGGCAGTGGGCAACCAAGTTGACTCGAGAAAGACGACCCTCCAGGATTTTCACTTGG
ACGAGGACAGGACTGTGAGAGTCCCATGATGTCAGACCTAAGGCCATTTACGATATGGCTTGGACTC
TGATCTCAACTGCAAGATTGCCAGCTGCCTTTGACAGGAAGCATGAGTATCATCTTCTCCTGCCCTG
ACGGTGACCCAGAATTGACCATGATAGAGGAGAGCCTCACCTCTGAGTTCGTTTCATGACATTGACCGG
AACTGAAGACTATCCAAGCTGTGCTGACTGTCCCAAGCTGAAGCTGAGCTATGAAGCGGACGTTACCAA
CTCTTTGACAGGACATGAAGCTACAGTCTTGTGTTGAGTCCCTGACTTCAGCAAGATTACCGGCAACCT
GTGAAGCTCACCAAGTGAACACAGGGCAGCTTTTGTGAGTGAATGAGGAGGGGGCAGGGACCAGCTCTA
ACCCAGACCTCCAGCCTGTCCGCTCACCTTCCCGCTCGACTATCACCTTAACCGACCGTTCATCTTTGT
TCTGAGGGACACGGACACGGGGCCCTCTTTCATAGGCAGAATCCTGGACCCAGCAGCACT**TAA**

ACGGTACGGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Chromatograms:	https://cdn.origene.com/chromatograms/ja1488_d08.zip
Restriction Sites:	Sgfl-Mlul
ACCN:	NM_177927
Insert Size:	1257 bp
OTI Disclaimer:	<p>Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.</p> <p>The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info</p>
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_177927.2 , NP_808788.1
RefSeq Size:	1434 bp
RefSeq ORF:	1257 bp
Locus ID:	287526
Cytogenetics:	10q24
Gene Summary:	mouse homolog is a proteinase inhibitor involved in vascularization [RGD, Feb 2006]