

Product datasheet for **RN215391**

Nfkbiz (NM_001107095) Rat Untagged Clone

Product data:

Product Type: Expression Plasmids

Tag: Tag Free

Symbol: Nfkbiz

Synonyms: Mail; RGD1310834

Mammalian Cell Selection: Neomycin

Vector: pCMV6-Entry (PS100001)

E. coli Selection: Kanamycin (25 ug/mL)

Fully Sequenced ORF: >RN215391 representing NM_001107095
Red=Cloning site Blue=ORF Orange=Stop codon

TTTGTGAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGC**

ATGATCGTGGACAAGCTGCTGGACGACAGCCGCGGAGAGGGGCTGCTGGACGCGCCGGCGACTGCG
GCCTCATGACCAGCCCGCTCAACCTCGCTTACTTCTACGGAGCCTCGCCTCCCTCCGCCCGCGCCGG
CGACACCGGTTACCTGTCCGCGGTTCCCTCCGCGCCCGGCTCGCCCGGCTCCGACTCTTCCGATTCTCC
TCCACTTCTCCGTGTCTCTGCGGGCCGTGGAGTCCCGTCCGAGAGGTGGCGCCCGCGCCGAGCGCC
CGCAAGTTGAGCCCATATGGGGGTGGCAGGCAGCAGAGAGGCCCTTTCAAGGTGTTTCGAGTGAAGAA
CTCGGTGAAAGAGCTGCTGCTGCACATCCGAAGTAATAAACAGAAGACTTCCGGCCAAGCGGTGGACGAG
TTTAAGAGTGCCTTTCAGGTGGCTGTGCTGCAATCAGCACCTCATTGTGCAAGACCTGGTGAACCTTG
GTGCCAGGTGAACACCACAGACTGCTGGGAAGGACACCTCTGCATGTCTGTGCAGAGAAGGGCCACTC
CCAGGTGCTCCAGGCAATCCAGAAGGGAGCAGTGAGGAGCAATCAGTTTGTGGATCTGGAGGCAACGAAC
TATGATGGCCTGACTCCCTCCACTGTGCGGTCTGGCTCACAATTCGGTGGTGACAGACCTGCAGAGAA
GCCAGCAGCCTATTGCGCCGAGGTGCAGGATCTTACTGAAGAATAAGAGCTTGGTTGACACCATTA
GTGCTGATTAGATGGAGCGCCGTGGAGGCAAGGATCGTAAAGTGGCCGACGGCCCTGCATTTG
GCAGCGGAGGAAGCAAACTCTGGAATCATTGCGCTCTTTTGGAGTTGCCGAGTTGCCTGTCTTTCGTGA
ATGCAAGGCTTACAATGGAACACCGCCCTCATGTTGCTGCCAGCCTGCAGTATCGGGTGACACAGTT
GGATGCAGTCCGTCTGTGATGAGGAAGGGAGCCGACCCAGTACTCGGAATGGAGAATGAACAGCCG
GTGCATTTGGTCCCTGATGGCCCTGTGGGAGAACAGATCCGACGTATCCTGAAGGGGAAGTCCATTACG
AGCGAGTACCACCATAC**TAG**

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



Restriction Sites:	Sgfl-Mlul
ACCN:	NM_001107095
Insert Size:	1140 bp
OTI Disclaimer:	Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery. Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>NM_001107095.1</u> , <u>NP_001100565.1</u>
RefSeq Size:	2990 bp
RefSeq ORF:	1140 bp
Locus ID:	304005
Cytogenetics:	11q12