

## Product datasheet for **RN215408**

### Edn1 (NM\_012548) Rat Untagged Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Edn1 (NM\_012548) Rat Untagged Clone  
**Tag:** Tag Free  
**Symbol:** Edn1  
**Synonyms:** Et1  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Cell Selection:** Neomycin  
**Fully Sequenced ORF:** >RN215408 representing NM\_012548  
**Red**=Cloning site **Blue**=ORF **Orange**=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**GCGATCGCC**

ATGGAATTATTTCCCGTGATCTTCTCTGCTGTTTGTGGCTTCCAAGGAGCTCCAGAAACAGCTGTCT  
TGGGAGCAGAGCTCAGCCCCGAGCTGAGAAGGAAGTGCAGAGCCCCCTCCAGCACATCTGGAGACC  
CCGCAGGTCCAAGCGTTGCTCCTGCTCCTTGTATGGACAAGGAGTGTGTCTACTTCTGCCACCTGGAC  
ATCATCTGGGTCAACACTCCCGAGCGCTCGTCCCGTATGGACTAGGAAGCCCTTCTAGGTCTAAGCGAT  
CCTTGAAAGACTTACTTCCCACAAAGACCACAGACCAAGGGAACAGATGCCAGTGTGCTCACAAAAAGA  
CAAGAAGTGTGGAATTTCTGCCAAGCAGACAAAGAACTCCGAGCCCAAAGTACCATGCAGAAAGGCGTA  
AAAGACTTCAAGAAGGGAAAACCTGTCCCAAGCTGGGAAAGAAGTGTATCTATCAGCAGCTGGTGGAGG  
GAAGAAAATAAGAAGTTGGAGGCCATCAGCAACAGCATCAAGACCTCCTTTCGAGTTGCAAAAGTTGAA  
AGCGGAAGTCTACAGAGACCAGAAGTTGATACACAACCGAGCACAT**TGA**

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Restriction Sites:** Sgfl-MluI  
**ACCN:** NM\_012548  
**Insert Size:** 609 bp  
**OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).



[View online »](#)

<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u><a href="#">NM_012548.2</a></u> , <u><a href="#">NP_036680.1</a></u>
<b>RefSeq Size:</b>	1390 bp
<b>RefSeq ORF:</b>	609 bp
<b>Locus ID:</b>	24323
<b>UniProt ID:</b>	<u><a href="#">P22388</a></u>
<b>Cytogenetics:</b>	17p13
<b>Gene Summary:</b>	activates ERK1/2 mediated signaling; may play a role in regulation of cardiac contractile force [RGD, Feb 2006]