

Product datasheet for **RN215413**

Sp6 (NM_001108833) Rat Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Sp6 (NM_001108833) Rat Untagged Clone
Tag:	Tag Free
Symbol:	Sp6
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>RN215413 representing NM_001108833 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGCTAACCGCTGTCTGTGGCTCTCTGGGCAGCCAGCACACCGACGCGCCTCACTCATCACCGCCGCGCC
TCGACCTGCAGCCTCTTCAGACATATCAGGGTCACACGAGCCCGGAGGCCGGGACTACCCCTCCCCGCT
GCAGCCTGGAGAGCTGCAGAGCCTCCGCTGGGCCCGGAGGTAGACTTCTCACAGGGCTATGAGTTGCCA
GGGGCCTCCTCGCGGGTAACCTGCGAGGACCTGGAAAGTGACAGTCCCTTGGCTCCGGGACCCCTTTCCA
AGCTCCTGCAGCCGACATGTACACCACTACGAATCGTGGTTCCGGCCGACGCACCCAGGCACGGAGGA
TGGCTCTTGGTGGGACCTTCATCCCGGCACCAAGCTGGATGGACCTCCCCATACTCAGGGAGCGCTGACC
TCACCTGGCCACCCGGGGCGCTTCAGCCTGCTTTGGGAGGATACGTGCGAGACCACCAAGCTCTGTGCTC
CGCCACCCACCCGACCCGACCACTCCTCCAGCCGCTGGTGGGCAGCACCTTTGGGGCCACCCGA
CGGGGCTAAGGTCTTGAAGCAGCAGCACCAGAGTCCCAAGGGCTGGATTCCAGTCTGGATGCCGCTGCC
CGACCCAAAGGCTCCCGGCGATCTGTGCCCCGAGCTCAGGGCAGACCGTGTGCTGCTGTCCCAACTGCC
TGGAGCGGAGCGACTCGGGCTCCGTGCGGGCCCGATGGGGCAAGAAGAAGCATTTCATAACTGCCA
CATCCAGGCTGCGGCAAGGCATACGCTAAGACGTCGATCTGAAGGCGCACCTGCGGTGGCACAGCGGC
GACCGACCTTCGTGTGAATTGGCTTTTCTGCGGCAAGCGCTTACGCGCTCCGACGAGCTGCAGCGCC
ATCTTCAGACCCACACCGGGACCAAGAAGTTCCCTGTGCGGTCTGCAGCCGTGTCTTCATGCGCAGCGA
CCACCTGGCCAAGCACATGAAAACCCACGAAGGCGCCAAAGAAGAGGCTGCGGCTGCGGGCCAGGGCGAG
GGCAAGGCCGGTGGCGTGGTGGAGTCACCCGGGGCAAAGGCAAACGCGAAGCCGAAGGCAGCTCGGCGT
CCTCCA**CTGA**

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI



[View online »](#)

ACCN:	NM_001108833
Insert Size:	1131 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001108833.1, NP_001102303.1</u>
RefSeq Size:	1995 bp
RefSeq ORF:	1131 bp
Locus ID:	363672
Cytogenetics:	10q31