

## Product datasheet for RN215461

### Crisp1 (NM\_001039393) Rat Untagged Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Crisp1 (NM\_001039393) Rat Untagged Clone  
**Tag:** Tag Free  
**Symbol:** Crisp1  
**Synonyms:** Crisp4  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Cell Selection:** Neomycin  
**Fully Sequenced ORF:** >RN215461 representing NM\_001039393  
 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCC**GCGATCGCC**

ATGGCAATGAAACTCATCTTGTTATTGTTTCTTGCTGCTACTCTTACTGTCTTTGTTCTGTTGTGACTC  
 TCAGACATCTCAAAGTGGATAGAGCTCTCTATAATCAACTGATCACGGAATCACAACTGAACCACAGGA  
 AGAAATTGTTGATACACACAATGCCTTCAGGAGAAACGTGTCCCGCCAGCCAGGAACATGCTGAAAATG  
 AGTTGGAGTTCAGCTGCTGCAGAAAATGCTAGAATTCTGGCGAGGTAAGTGCAGACAAGTGCAGACGTGACT  
 CACTTGAAAGGAGACTACCAAATACATTTTGTGGAGAAAATATGCATATGGAAAACATCCCTCTTCTCTG  
 GTCAAATGTAATTGAAATCTGGTACAATGAGTCGAAAATATTTCAAATACGGAGAAATGGCCATCCACAGAT  
 GATGACATTGAACTTATCACTACACCCAGATGGTCTGGGCTCTTCTACCTCATTGGCTGTGATGTTG  
 CATCATGCCCGAGGCAAAAGGCAGCTACGTATCTCTATGTGTGCTACTATTGCCATGAGGGGAATAGTCA  
 AGACACATTAACATGCCTTATAAGGAGGGGCCCATGTCAAGACTGTCCAAACAACGTGAAGATGGA  
 CTGTGTACCAACCCTGCCTTTATTACGATGAATATAACAACGTGATAAGCAAGTAAAACCTCTGGAT  
 GCTCACATCCAGCAGTCCGGCCATTCTGTAAGGCTTTTGTCTGTGTA AAACTGAAATAAAAT**TAA**

**ACGCGT**ACGCGGCCGCTCGAGCAGAAAACATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
 ACAAGGATGACGACGATAAGGTTTAA

**Restriction Sites:** Sgfl-Mlul  
**ACCN:** NM\_001039393  
**Insert Size:** 765 bp



[View online »](#)

<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u><a href="#">NM_001039393.2</a></u> , <u><a href="#">NP_001034482.2</a></u>
<b>RefSeq Size:</b>	1394 bp
<b>RefSeq ORF:</b>	765 bp
<b>Locus ID:</b>	654517
<b>Cytogenetics:</b>	9q13