

Product datasheet for RN215716

Uqcr10 (NM_001170465) Rat Untagged Clone

Product data:

Product Type: Expression Plasmids

Product Name: Uqcr10 (NM_001170465) Rat Untagged Clone

Tag: Tag Free

Symbol: Uqcr10

Synonyms: Ucrc

Mammalian Cell Selection: Neomycin

Vector: pCMV6-Entry (PS100001)

E. coli Selection: Kanamycin (25 ug/mL)

Fully Sequenced ORF: >RN215716 representing NM_001170465
 Red=Cloning site Blue=ORF Orange=Stop codon

TTTGTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCCGCGATCGCC

ATGGCGTCGCCGACGGTCACTTCGCGTTTGTACTCCTTGCTGTTCCGCAGGACTTCTACCTTCGCCCTCA
 CCATCGCAGTGAGCGCCCTGTTCTTCGAGCGCGCCTTCGATCAGGGTGCAGACGCGGTCTACGATTACAT
 CAACGAGGGGAACTGTGGAACACATAAAGCACAAAGTACGAGAAACAGGAGTAA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI

ACCN: NM_001170465

Insert Size: 195 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: Clone contains native stop codon, and expresses the complete ORF without any c-terminal tag.



[View online »](#)

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001170465.1, NP_001163936.1</u>
RefSeq Size:	455 bp
RefSeq ORF:	195 bp
Locus ID:	685322
Cytogenetics:	14q21