

Product datasheet for **RN215798**

Dpy30 (NM_001170545) Rat Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Dpy30 (NM_001170545) Rat Untagged Clone
Tag: Tag Free
Symbol: Dpy30
Synonyms: Aip1
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >RN215798 representing NM_001170545
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**GCGATCGCC**

ATGGAGTCAGAGCAGATGCTGGAGGGACAGACACAGGTTGCAGAAAACCTCACTCTGAGTACGGGCTCA
CAGACAGCGTTGAGAGAATAGTTGAAAATGAGAAGATTAATGCAGAAAAGTCATCAAACAGAAGGTGGA
TCTACAGTCCTTGCCACCCGCTACCTGGACCAGACAGTTGTGCCTATCTTATTACAGGGACTTGCT
GTGCTCGCAAAGGAAAGACCACCAAATCCCATCGAGTTTCTAGCATCTTATCTTTAAAAAATAAGGCAC
AGTTTGAAGATCGAAAT**TGA**

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: Sgfl-MluI
ACCN: NM_001170545
Insert Size: 300 bp
OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation: Clone contains native stop codon, and expresses the complete ORF without any c-terminal tag.



[View online »](#)

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001170545.1, NP_001164016.1</u>
RefSeq Size:	680 bp
RefSeq ORF:	300 bp
Locus ID:	286897
Cytogenetics:	6q13
Gene Summary:	<p>As part of the MLL1/MLL complex, involved in the methylation of histone H3 at 'Lys-4', particularly trimethylation. Histone H3 'Lys-4' methylation represents a specific tag for epigenetic transcriptional activation. May play some role in histone H3 acetylation. In embryonic stem cells, may play a crucial role in retinoic acid-induced differentiation along the neural lineage, regulating gene induction and H3 'Lys-4' methylation at key developmental loci. May also play an indirect or direct role in endosomal transport (By similarity). [UniProtKB/Swiss-Prot Function]</p> <p>Transcript Variant: This variant (2) differs in the 5' UTR, compared to variant 1. Variants 1 and 2 encode the same protein.</p>