

## Product datasheet for RN215824

### Pate2 (NM\_001267781) Rat Untagged Clone

#### Product data:

|                           |                                                                                    |
|---------------------------|------------------------------------------------------------------------------------|
| Product Type:             | Expression Plasmids                                                                |
| Product Name:             | Pate2 (NM_001267781) Rat Untagged Clone                                            |
| Tag:                      | Tag Free                                                                           |
| Symbol:                   | Pate2                                                                              |
| Synonyms:                 | Pate2                                                                              |
| Mammalian Cell Selection: | Neomycin                                                                           |
| Vector:                   | pCMV6-Entry (PS100001)                                                             |
| E. coli Selection:        | Kanamycin (25 ug/mL)                                                               |
| Fully Sequenced ORF:      | >RN215824 representing NM_001267781<br>Red=Cloning site Blue=ORF Orange=Stop codon |

TTTGTAAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCC**CGATCGCC**

ATGCTGGTGATAGTCTGTTTGTCTGCCAGTATTGGGGTGTCTTAATGAACCTGAAAAAGATCTTGAA  
 CAATGTGTTATAAATGAAGAAATATCATCTTGGGTTATGCTACGGACTCATGAGATCTGCACACTGAA  
 GCATAGACAGTCCTGTGCTGCTGAGAACTTTACATACTACGAACAGAGGGCAGAGCATGTATCATTAT  
 TCAAGACTGTATGTATGACCAACTGTGAGGACATCAACTTCCTGAGTTTGAAGGAGGACAGAGCTAA  
 TTTGTTGCAAGCACAGTAGCTATTGCAACCTCCAATGGGACT**C**TA**G**

**ACGCGT**ACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
 ACAAGGATGACGACGATAAGGTTTAA

|                    |              |
|--------------------|--------------|
| Restriction Sites: | Sgfl-MluI    |
| ACCN:              | NM_001267781 |
| Insert Size:       | 327 bp       |

**OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

**OTI Annotation:** Clone contains native stop codon, and expresses the complete ORF without any c-terminal tag.


[View online »](#)

|                               |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |
|-------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Components:</b>            | The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).                                                                                                                                                                                                                                                                                                                              |
| <b>Reconstitution Method:</b> | <ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol> |
| <b>RefSeq:</b>                | <u>NM_001267781.1</u> , <u>NP_001254710.1</u>                                                                                                                                                                                                                                                                                                                                                                                                                                                                     |
| <b>RefSeq Size:</b>           | 327 bp                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                            |
| <b>RefSeq ORF:</b>            | 327 bp                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                            |
| <b>Locus ID:</b>              | 100361465                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
| <b>Cytogenetics:</b>          | 8q21                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                              |