

Product datasheet for RN215832

Eloc (NM_001270563) Rat Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Eloc (NM_001270563) Rat Untagged Clone
Tag:	Tag Free
Symbol:	Eloc
Synonyms:	Tceb1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>RN215832 representing NM_001270563 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGATGGAGAGGAGAAAACCTATGGTGGCTGTGAAGGCCCTGATGCCATGTATGTGAAATTAATATCTT
CTGATGGTCATGAATTTATTGTAAAAAGAGAACATGCATTAACATCAGGAACAATAAGGCCATGTTGAG
TGGTCCAGGTCAGTTTGGGAGAATGAAACCAATGAGGTCAACTTAGAGAGATCCCTTCACATGTGCTA
TCGAAAGTGTGCATGTATTTACCTACAAGGTCCGCTATACTAACAGCTCCACTGAAATTCCTGAATTCC
CAATTGCACCTGAAATTGCACTGGAAGTCTGATGGCCGCGAATTCCTAGATTGT**TAA**

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites:	Sgfl-Mlul
ACCN:	NM_001270563
Insert Size:	339 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: Clone contains native stop codon, and expresses the complete ORF without any c-terminal tag.



[View online »](#)

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001270563.1</u> , <u>NP_001257492.1</u>
RefSeq Size:	2014 bp
RefSeq ORF:	339 bp
Locus ID:	64525
UniProt ID:	<u>P83941</u>
Cytogenetics:	5q11
Gene Summary:	<p>15 kDa subunit of the transcription factor B (SIII) complex, also known as elongin C; releases RNA polymerase II from transient pausing at arresting sites [RGD, Feb 2006]</p> <p>Transcript Variant: This variant (4) differs in the 5' UTR compared to variant 1. All four variants encode the same protein. Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.</p>