

Product datasheet for **RN215875**

Fau (NM_001160231) Rat Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Fau (NM_001160231) Rat Untagged Clone
Tag: Tag Free
Symbol: Fau
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >RN215875 representing NM_001160231
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGCAGCTCTTTGTCCGCGCCAGGAACACTACACCCCTCGAGGTGACCGCCAGGAGACGGTCGCCCAGA
TCAAAGCTCATGTGGCCCTCACTGGAAGGCATTGCCCCGAAGATCAAGTCGTGCTTCTGGCAGGCTCGCC
GCTGGAGGATGAGGCCACCCTAGGCCAGTGTGGTGTAGAGGCCCTGACCACTCTGGAGGTAGCTGGCCGC
ATGCTGGGAGGTAAGTTTCATGGTCTCTGGCTCGTGCTGGGAAAGTGAGAGGTCAGACTCCAAGGTGG
CCAACAGGAAAAGAAGAAGAAGACTGGTCCGGCCAAGAGGCGAATGCAGTACAACCGACGCTTTGT
CAATGTTGTGCCACCTTTGGCAAGAAGAAGGCCCAATGCTAATTCT**TAA**

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI
ACCN: NM_001160231
Insert Size: 402 bp
OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation: Clone contains native stop codon, and expresses the complete ORF without any c-terminal tag.



[View online »](#)

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001160231.1</u> , <u>NP_001153703.1</u>
RefSeq Size:	545 bp
RefSeq ORF:	402 bp
Locus ID:	29752
Cytogenetics:	1q43
Gene Summary:	hybrid protein; encodes a ubiquitin-like protein and ribosomal S30 protein [RGD, Feb 2006] Transcript Variant: This variant (2) differs in the 5' UTR, compared to variant 1. Variants 1-3 encode the same protein.