

Product datasheet for RN215897

Lyzl4 (NM_001246183) Rat Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Lyzl4 (NM_001246183) Rat Untagged Clone
Tag: Tag Free
Symbol: Lyzl4
Synonyms: Lyza; RGD1308401
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Fully Sequenced ORF: >RN215897 representing NM_001246183
 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAAACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGCAGCTGTACCTGGTCTTCTCCTCATCAGCTACTTGCTGACTCCGATTGGTGCCTCCATTTGGGGC
 GCTGCGTGGTGGCTAAGAACTCTATGACGGGGGCTTAAATTATTTGAGGGCTACAGCCTTGAAACTG
 GGTGTGCCTGGCTATTTTGAAGCAAGTTCAACCCCTCGGCTGTCTATGAGAACTACGAGATGGCTCC
 ACCGGCTTCGGCTCTTTCAGATCCGCGACAATGAATGGTGTGACCATGGCAAGAACCTCTGCAGTGTGT
 CCTGCACGGCTTTACTGAATCCCACTTAAAGGATACAATCGAGTGTGCCAAGAAAATCGTGAAAGGAAA
 ACAAGGGATGGGAGCATGGCCCGTCTGGTCCAGGAAGTCCAGCTGTCCGACATTCTGGATCGGTGGCTG
 GATGGCTGTGAGCTG**TAG**

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI
ACCN: NM_001246183
Insert Size: 438 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).


[View online »](#)

OTI Annotation:	Clone contains native stop codon, and expresses the complete ORF without any c-terminal tag.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001246183.1</u> , <u>NP_001233112.1</u>
RefSeq Size:	580 bp
RefSeq ORF:	438 bp
Locus ID:	363168
UniProt ID:	<u>D4ABW7</u>
Cytogenetics:	8q32
Gene Summary:	May be involved in fertilization (By similarity). Has no detectable bacteriolytic and lysozyme activities in vitro (PubMed:22110709).[UniProtKB/Swiss-Prot Function]