

Product datasheet for **RN215901**

Wfdc6a (NM_001246283) Rat Untagged Clone

Product data:

Product Type: Expression Plasmids

Product Name: Wfdc6a (NM_001246283) Rat Untagged Clone

Tag: Tag Free

Symbol: Wfdc6a

Mammalian Cell Selection: Neomycin

Vector: pCMV6-Entry (PS100001)

E. coli Selection: Kanamycin (25 ug/mL)

Fully Sequenced ORF: >RN215901 representing NM_001246283
 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGCCTGCCAACAGGTTACAGCTTCCAAAATGAGACACTGGGGGCTTCTGCCATTCTGGTCCCCCTCA
 TCTTCCTGTGGAGCATCCAGAACTGAGCTGACAGAGGGTTATTTTCAGATCATGTCCGAGAGTAAG
 AGTGAAGTGTGAGGTTGAAGAAAGAAATGAGTGTACGAGGCACCGACAGTGTCCGGACAAGAAGAGATGC
 TGCTTCTTCAGCTGTGGGAAGAAGTGCATGGACCTCCGGAAGACATATGCAGTCTGCCACAGGATGCTG
 GCCCTGCCTGGCTACCTCCACGCTGGTGGTATAATGAAGATACTGGCCTCTGTACCCAATTCATCTA
 TGGTGGTTGCCAAGGGAACCCCAACAATTCCAATCTGAAGGTATCTGCACAGTTGTCTGCAAAAAGAAA
 CAAATGTCTTCTGGAT**ATGA**

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI

ACCN: NM_001246283

Insert Size: 441 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).



[View online »](#)

OTI Annotation:	Clone contains native stop codon, and expresses the complete ORF without any c-terminal tag.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001246283.2, NP_001233212.2</u>
RefSeq Size:	749 bp
RefSeq ORF:	441 bp
Locus ID:	685153
Cytogenetics:	3q42