

## Product datasheet for **RN215906**

### **Atoh7 (NM\_001170482) Rat Untagged Clone**

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Atoh7 (NM\_001170482) Rat Untagged Clone  
**Tag:** Tag Free  
**Symbol:** Atoh7  
**Mammalian Cell Selection:** Neomycin  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Fully Sequenced ORF:** >RN215906 representing NM\_001170482  
**Red**=Cloning site **Blue**=ORF **Orange**=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGAAGTCCGCCTGCAAACCCACGGTCTCCGGCGGGAGCACGCGGCGGCCCCCGTGC GCGGGCGCAG  
CCGAGCGCGCGGTCTCTGCGCGGGGCCGGCGGCTGAAAGCGCAGCGCGCAGGCGCCTGGCGGCCAA  
CGCGCGGAGCGGCGCCGATGCAAGGGCTGAACACCGCGTTCGACCGGCTGCGCAGGGTGGTGCCGCAA  
TGGGGCCAGGACAAGAAGCTGTCCAAGTACGAGACACTGCAGATGGCGCTCAGCTACATCGTGGCGCTCA  
CCCGATCCTAGCCGAAGCCGAGCGGGACTGGGTGGGGCTGCGCTGCGAGCAGCGGGCCGCAATCACCC  
CTACCTCCATTCCCGGTGCTAGGCTCCAGGTAGACCCTGAGCCCTACGGGCAGAAGTTCTTCGGTTTC  
CAGCCAGAGTCTTCCCATGGCCAGCTAA

**ACGCGT**ACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Restriction Sites:** SgfI-MluI  
**ACCN:** NM\_001170482  
**Insert Size:** 450 bp



[View online »](#)

**OTI Disclaimer:** Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at [custsupport@origene.com](mailto:custsupport@origene.com) or by calling 301.340.3188 option 3 for pricing and delivery.

Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

**OTI Annotation:** Clone contains native stop codon, and expresses the complete ORF without any c-terminal tag.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_001170482.1](#), [NP\\_001163953.1](#)

**RefSeq Size:** 450 bp

**RefSeq ORF:** 450 bp

**Locus ID:** 365564

**Cytogenetics:** 20p11