

Product datasheet for RN215945

LOC100363531 (NM_001177826) Rat Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: LOC100363531 (NM_001177826) Rat Untagged Clone
Tag: Tag Free
Symbol: LOC100363531
Synonyms: LRRGT00182
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Fully Sequenced ORF: >RN215945 representing NM_001177826
 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAAACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGGAGGCTGGTGACATAGCTCAGCAGCAGAGCTTTTGTGTAGGATATACAATCCAGCACAACACACA
 CGCACACACGCACACACACACACAGAAATAGAAAATACAGGATCCACTTCAGACTAAAGTTGA
 AACAGGACTAATGGCATCCATCATTTAACAGTACTTCAGAGTCCCTCAATCTTAAAGATGCTAATGCT
 TCTAAGGTCTCTGAAGAAAGAGAATGCAGAACTAGGCTAGGAACACAAGGACAGGATGATAAGCCACTG
 TCACAGTACCACTGTCATTGCTATCTTCCAGTACTGGAGACCACGCAGGGAACCGTGCACCTAGAAT
 TAAAACTGCAGGCTTTAGGGCAGAGGAAAAAGGCATGAGAGCACGAGGCAGATGGATCCCGGGTGTGT
 AAAGGGCTCGCCACATTATCTGCTCACCTTGGAGACGACTTAGTGCGTTTTGGCCAAGGGCAACATCG
 CGGAGACAAGGATGGAAGTGC**GTGA**

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI
ACCN: NM_001177826
Insert Size: 516 bp
OTI Disclaimer:

Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).


[View online »](#)

OTI Annotation:	Clone contains native stop codon, and expresses the complete ORF without any c-terminal tag.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001177826.1</u> , <u>NP_001171297.1</u>
RefSeq Size:	516 bp
RefSeq ORF:	516 bp
Locus ID:	100363531
Cytogenetics:	Xq22