

Product datasheet for RN215956

Ascl3 (NM_001271234) Rat Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Ascl3 (NM_001271234) Rat Untagged Clone
Tag:	Tag Free
Symbol:	Ascl3
Synonyms:	Sgn1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>RN215956 representing NM_001271234 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGACACCAGAAGTTATCCCAGCCCTCCAGACAGGCTCTCGGTCTTCGCAGAGTCTGCCCACTTGCCAC
TGTCCAGGCCCTTCTACCTGGACCCCATGGTTACTGTCCACCTATGCCCCGACGCCCGGTGCCGGCCCC
TTACACGGACGAGCTGCCTCTGCTGCCCTTTCCAGTGACACCCTGATCATGGACAATTACGGGGACCCA
TACCCCTTCCCCTTCCCTATGCCTTACACCAACTACAGGGCTGTGACTATACCTATGGCCGGCCTTCA
TCCGCAAGAGGAACGAGCGAGAGAGGCAGCGAGTCAAGTGTGTGAACGAAGGCTATGCCCGCTGCGACG
GCACCTGCCAGAGGACTACCTGGAGAAACGACTCAGCAAAGTGGAGACCCTCAGAGCTGCCATCAAGTAT
ATCAGCTACCTGCAGTCTCTTGTACCCGGATGAGACTGAGACCAAGAAAACCCTCGAACAACCAACT
GTGGCTCCCTGGACCCGGCTCTGAGAGTCATT**TGA**

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGTTTAA

Restriction Sites:	Sgfl-MluI
ACCN:	NM_001271234
Insert Size:	525 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).



[View online »](#)

OTI Annotation:	Clone contains native stop codon, and expresses the complete ORF without any c-terminal tag.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001271234.1</u> , <u>NP_001258163.1</u>
RefSeq Size:	689 bp
RefSeq ORF:	525 bp
Locus ID:	246301
Cytogenetics:	1q33
Gene Summary:	mouse homolog represses transcription and may play a role in salivary gland development [RGD, Feb 2006]