

Product datasheet for RN215959

Bles03 (NM_001243534) Rat Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Bles03 (NM_001243534) Rat Untagged Clone
Tag: Tag Free
Symbol: Bles03
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Fully Sequenced ORF: >RN215959 representing NM_001243534

Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGGCGGCGCGCGCGCGCGCGTGGCAGGGGCGGGGCGCGCGGGGGTGGCGGCGCAGACCCGGGCAGG
 AACGGAGCCGGGCCGAAGTTGGGTTGGCGCGGAGCGAGTGAAGGCCGAGTTCCTCACAACAGGATG
 GAACCAAAATGAAGAGCTGGAAGAGGAGGACTCTCCAGGTGGTCGTGAGGATGGCTTCACTGCTGAGCACC
 TGGCTGCAGAGGCCATGGCAGCTGACATGGACCCCTGGCTGGTATTTGATGCCCGTACTACACCTGCCAC
 AGAGCTGGATGCCTGGTTGGCCAAGTACCCACCATCTCAAGTTACTCGCTATGGGGACCCAGGTTACCC
 AACTCTGAACCTGTGGGCTGGATTGCAGCCTATGGGCAGGGTTACACCCCAATTACAGGGGATGTACAAG
 GGCTGCAGGCAGCCTGGGAGGCTCTGCAGACCAGTGGGCGACCCATCACACCAGGTACCCTGCGCCAGCT
 GGCCATCACCCACCATGTGCTCTCCGCAAGTGGC**TGA**

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI
ACCN: NM_001243534
Insert Size: 528 bp
OTI Disclaimer:

Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).



[View online »](#)

OTI Annotation:	Clone contains native stop codon, and expresses the complete ORF without any c-terminal tag.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001243534.1</u> , <u>NP_001230463.1</u>
RefSeq Size:	1613 bp
RefSeq ORF:	528 bp
Locus ID:	266609
UniProt ID:	<u>Q566Q8</u>
Cytogenetics:	1q43