

Product datasheet for **RN215972**

Il17b (NM_053789) Rat Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Il17b (NM_053789) Rat Untagged Clone
Tag: Tag Free
Symbol: Il17b
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Fully Sequenced ORF: >RN215972 representing NM_053789
Red=Cloning site **Blue**=ORF **Orange**=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGACTGGCCACACAGCCTGCTCTTCCTTCTTGCCATCTCCATCTTCCTGGTGCCAAGCCAGCCCCGGA
ACACCAAAGGCCAAAAGGAAGGGGCAAGGGAGGCCCTGGCCCTTGGCCCTGGCCCTACCCAGGTGCCGCT
GGACCTGGTGTCTCGAGTAAAGCCCTACGCCGAATGGAAGAGTACGAGAGGAACCTTGGCGAGATGGTG
GCCAGCTGAGGAACAGCTCTGAGCCAGCCAAGAAGAAGTGTGAAGTCACTTGCAGCTGTGTTATCCA
ACAAGAGAAGCCTGTCTCCCTGGGGCTACAGCATCAATCAGACCCAGCCGCATCCCGGAGGACTTGCC
TGAGGCGCGGTGCCTATGTTGGGTTGCGTGAACCCCTTACCATGCAGGAGGACCGTAGCATGGTGAGC
GTGCCAGTGTTAGCCAGGTGCCAGTGCCGCCGCCCTCTGTCCGCAACCTCCTCGGCCGGCCCTGCC
GCCAGCGTGTTGTCATGGAGACCATCGCTGTGGGTTGCACCTGCATCTT**CTGA**

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI
ACCN: NM_053789
Insert Size: 543 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).



[View online »](#)

OTI Annotation:	Clone contains native stop codon, and expresses the complete ORF without any c-terminal tag.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_053789.1</u> , <u>NP_446241.1</u>
RefSeq Size:	688 bp
RefSeq ORF:	543 bp
Locus ID:	116472
Cytogenetics:	18q12.1
Gene Summary:	human homolog is a T cell-derived cytokine that may be involved in the proinflammatory response [RGD, Feb 2006]