

Product datasheet for RN215973

Spcs3 (NM_001191073) Rat Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Spcs3 (NM_001191073) Rat Untagged Clone
Tag:	Tag Free
Symbol:	Spcs3
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin
Fully Sequenced ORF:	>RN215973 representing NM_001191073 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGAACACCGTGTCTCGCGCGGAATTCCTGTTTCGCTTCTCGCTGAGTGTGATGGCGCGCTCACCT
TCGGCTGTTTCATCACCACCGCTTCAAAGACAGGAGCGTCCCGGTGCGGCTGCACGTCTCGCGGATCAT
GCTAAAAAATGTAGAAGACTTCACTGGCCCTAGAGAAAGAAGTGACCTTGGATTCACATTTGATATA
ACAGCTGATCTAGAGAATATATTTGATTGGAATGTTAAGCAATTATTTCTTTATTTATCAGCAGAATATT
CAACAAAAAATAATGCTCTGAACCAAGTTGTCCTTTGGGACAAGATTGTTTTGAGAGGTGATAATCCAAA
GCTGCTTTTGAAGATATGAAGCAAAAGTATTTCTTTGACGATGGAAATGGCCTCAAGGAAACAGG
AATGCTACTTTGACGCTCTCTGGAATGTTGTACCAAACGCTGGAATTCTACCTTTGTGACAGGATCGG
GACATATATCTGTCCATTTCCAGATACATACGAAATCACCAAGAGTTACT**AA**

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites:	Sgfl-MluI
ACCN:	NM_001191073
Insert Size:	543 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	Clone contains native stop codon, and expresses the complete ORF without any c-terminal tag.



[View online »](#)

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001191073.1, NP_001178002.1</u>
RefSeq Size:	763 bp
RefSeq ORF:	543 bp
Locus ID:	680782
Cytogenetics:	16p11