

Product datasheet for **RN215984**

RGD1563365 (NM_001271340) Rat Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	RGD1563365 (NM_001271340) Rat Untagged Clone
Tag:	Tag Free
Symbol:	RGD1563365
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>RN215984 representing NM_001271340

Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGCGCTCCGTGGTGAGCACGCGAGTGGTCTGGGGTTTCTTCGTCTTCTGTGGACTGCTCGCGGATCA
CGGTGTCCAGGGAGCTACTGACGGCGGGGAGCGAGGGTTCGCGAGGTGTCTGGGACCAGCTGCTTATCAA
CTCTAAGCCTCATTCCAGAAAGGCCTCTACCCTTCAAACAGTTCGGGTACAAAGGAGCCCTTTGTTAGAC
CAGGTTTCAGGCCTTCTCCACAAATGGCCAGGCAATGAGAAGCTCGGAAGGAAATGGCTGTTGCTC
CAGCTGGTCATTCAATATTGAAAACATCGACGAGACGTCCGGAACATCATAACAGTGGACGTGGCCTT
ATTTGAGATGAATCAGTCAGAAGATTGAAAGAGGACGATAATCCAGAAGAGAGCTCACAGGACAGTTCT
GGGGACAGTTCAGAGTCTGAAGACGACGACTGCATTCCCTCTGAAGTCACCATCGAGAACATTAAGCTTC
CGAACACAGAAGGTGGCAAAGGCAAGATTCAAGTTTGGACAGTCCGGCGAGTAATAAAACAATCAGTC
AAAT**TAA**

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites:	Sgfl-MluI
ACCN:	NM_001271340
Insert Size:	567 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).



[View online »](#)

OTI Annotation:	Clone contains native stop codon, and expresses the complete ORF without any c-terminal tag.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001271340.2</u> , <u>NP_001258269.1</u>
RefSeq Size:	1189 bp
RefSeq ORF:	567 bp
Locus ID:	299700
Cytogenetics:	7q13