

Product datasheet for **RN215990**

Park7 (NM_001277251) Rat Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Park7 (NM_001277251) Rat Untagged Clone
Tag:	Tag Free
Symbol:	Park7
Synonyms:	CAP1; DJ-1; Dj1; SP22
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin
Fully Sequenced ORF:	>RN215990 representing NM_001277251 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**GCGATCGCC**

ATGGCATCCAAAAGAGCTCTGGTCATCCTAGCCAAAGGAGCAGAGGAGATGGAGACAGTGATTCTGTGG
ACATCATGCGCGAGCTGGGATTAAGTCACCGTTGCAGGCTTGGCTGGGAAGGACCCCGTCAGTGTAG
CCGTGATGTAGTATTTGTCCGGATACCACTCTGGAAGAAGCAAAAACACAGGGACCATACGATGTGGTT
GTCTTCCAGGAGGAAATCTGGGTGCACAGAATTATCTGAGTCGGCTTTGGTGAAGGAGATCCTCAAGG
AGCAGGAGAACAGGAAGGGCCTCATAGCTGCCATCTGTGCGGGTCTACGGCCCTGCTGGCTCACGAAGT
AGGCTTTGGATGCAAGGTTACATCGCACCCATTGGCTAAGGACAAAATGATGAACGGCAGTCACTACAGC
TACTCAGAGAGCCGTGTGGAGAAGGACGGCCTCATCTCACCAGCCGTGGGCCTGGGACCAGCTTCGAGT
TTGCGCTGGCCATTGTGGAGGCACTCAGTGGAAGGACATGGCTAACCAAGTGAAGGCCCCGCTTGTCT
CAAAGC**TAG**

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites:	Sgfl-MluI
ACCN:	NM_001277251
Insert Size:	570 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).



[View online »](#)

OTI Annotation:	Clone contains native stop codon, and expresses the complete ORF without any c-terminal tag.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001277251.1</u> , <u>NP_001264180.1</u>
RefSeq Size:	931 bp
RefSeq ORF:	570 bp
Locus ID:	117287
UniProt ID:	<u>O88767</u>
Cytogenetics:	5q36
Gene Summary:	displays a specific localization pattern to the anterior-ventral surface of the equatorial segment; may play a role in sperm-egg interactions and fertilization [RGD, Feb 2006] Transcript Variant: This variant (3) differs in the 5' UTR and includes an alternate penultimate exon, compared to variant 1, resulting in a novel 3' coding region and 3' UTR. The encoded isoform (2) is shorter, compared to isoform 1. Variants 2, 3, 4, 5, and 6 encode the same isoform. Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.