

Product datasheet for RN216028

Msc (NM_001191754) Rat Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Msc (NM_001191754) Rat Untagged Clone
Tag: Tag Free
Symbol: Msc
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >RN216028 representing NM_001191754
 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGTCCACCGGCTCGGTGAGTGACCCCGAAGACACGGAGATGCGGGGGCTGCAGCGGGTCTACCCGGCCC
 CCGCCTCCAAGAGGCCACCCCTGCTCCGCCCGGAGCGCAGTTACGCCCTGCCAGCGACAACCTTCGGC
 GGAGGAGGAGACCCCGACGGCGAAGAGGAGCCGGCACCGTGGGAGCCCGGGAGGCTGCAAGAGAAAAG
 CGGCCCGTGGGGCTGACGCCAGCGCGCAGGTGGCGGTGCGGGCAGCACCGGGAAGAAGCGCTCCCGC
 CCAAGGGCTCGGCCCGGAGTGCAAGCAGTCACAGCGGAACGCGCCAACGCCCGCAACGCGCCCGGAT
 GCGCGTGTGAGCAAAGCTTTCTCCAGACTCAAGACCAGCCTGCCTTGGGTGCCGCCGACACCAAGCTT
 TCCAAACTGGACACGCTGCGTTTGGCTTCCAGCTACATCGCGCACCTGCGCCAGCTGCTGCAGGAGGACC
 GCTACGAGGATAGCTATGTGCACCCCGTGAACCTGACGTGGCCATTCTGTGGTCTCTGGACGCCAGACTC
 TGATACCAAAGATGTTTCTGCAGCCAACAGGCTTTGTGGAACCTCAGCT**TAG**

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI
ACCN: NM_001191754
Insert Size: 612 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).



[View online »](#)

| | |
|-------------------------------|---|
| OTI Annotation: | Clone contains native stop codon, and expresses the complete ORF without any c-terminal tag. |
| Components: | The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water). |
| Reconstitution Method: | <ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C. |
| RefSeq: | <u>NM_001191754.1, NP_001178683.1</u> |
| RefSeq Size: | 612 bp |
| RefSeq ORF: | 612 bp |
| Locus ID: | 312897 |
| Cytogenetics: | 5q11 |