

Product datasheet for RN216045

LOC100363112 (NM_001177822) Rat Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: LOC100363112 (NM_001177822) Rat Untagged Clone
Tag: Tag Free
Symbol: LOC100363112
Synonyms: Aa2-296
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Fully Sequenced ORF: >RN216045 representing NM_001177822
 Red=Cloning site Blue=ORF Orange=Stop codon

TTTGTAAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGCCACAGCCTCAAAGAAGCCCGAGTTACCTTACGGCGACACAAATAAGAACTCAATCTGCTCACC
 CCAGTGTAGGACAGTGTGAAGAGGACCGTTGTATTTGACGGTGTAGGAAATGGAATTCACAGGGCTT
 TCCGTAATCACAATGGGCACCAAGACCCACTGCTGCACAAGTAGACGCCAACGGAAGTGAAGAGGAAG
 TATGTAGAAAGACCCTGTCTCAAAGAGGAGGAGGAGGAGGAAGAGGAGGAAGAAGAGGAGGAGAAGA
 AAGAAGGAGAAGAAGAAGAAGAAGAAGAAGAAGAAGAAGAAGAAGAAGAAGAAGAAGAAGAAGA
 AGAAGAAGAAGAAGAAGAAGAAGAAGAAGACACAAAGACAGAAATAAATAAGACACTGGACCACCAACCA
 GGCAGCACCACCGCCTCTGTGAGGCCCGGACACATACATAGCAGAGGACCGCCAGGTCTGGGTACAGT
 CAGAGAAGATGCACCTAACCCTCAGGAGACTGGAGGCCCGAGGAGTGGGGAGGGGTGGTGGTGGGAAT
 CCTCATGGAGACGGGTGGGGTGGTATGGGATGTGGGCAGAGGGTAGATGAAGATAAAATGA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI
ACCN: NM_001177822
Insert Size: 624 bp


[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	Clone contains native stop codon, and expresses the complete ORF without any c-terminal tag.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001177822.1, NP_001171293.1</u>
RefSeq Size:	670 bp
RefSeq ORF:	624 bp
Locus ID:	100363112
Cytogenetics:	3q36