

Product datasheet for **RN216051**

Igll1 (NM_001190341) Rat Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Igll1 (NM_001190341) Rat Untagged Clone
Tag: Tag Free
Symbol: Igll1
Synonyms: Igll1
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >RN216051 representing NM_001190341
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGAAGCTCAGGGCAGGACAGACACTGGGCACTATCCCCAGGCAATGTGAAATTCTCCTTCTGTTGCTGC
TGTTGGGCTGGTGGATGGTGTCCACCATATACTTTCCCCAAGCTCAGCACAAAGGGGCAGAGCTGTGGG
CCCTGGAGCCTCAGTGGGAAGCAGCAGGTCTAGCCTGTGGACCCTCCAGGCAGGTTCTGTTCCAGATC
ATCCACGGGGAGCAGGTCCAGGTGCTGGCCCCATAGGCTTCCATCCAAGTCCCAGTTATGGTAGTCT
TTGGTAGTGGGACCCAGCTCACAATCCTAGGGCAGCCCAAGTCTGACCCCTTGGTCACTCTGTTCTGCC
TTCCTTAAAGAATCTCCAGGTTAAGAAGGCGACACTAGTGTGTCTGGTGAGCGAATTCTACCCAGGTA
TTGGTGGTGGACTGGAAGGTAGATGGGATCCCTGTCACTCAGGGTGTGGAGACAACCAACCTCCAAAC
AGACCAACAACAAGTACGTGGCCAGCAGCTACCTGACACTGATGTCTGACCAATGGATGCCTCACAGTAG
ATACACCTGCCAGGTCACTCATGAAGGAAACACTGTGGAGAAGAGTGTGTACCTGCTGAATGTTCT**TAG**

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI
ACCN: NM_001190341
Insert Size: 630 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).



[View online »](#)

OTI Annotation:	Clone contains native stop codon, and expresses the complete ORF without any c-terminal tag.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001190341.1</u> , <u>NP_001177270.1</u>
RefSeq Size:	630 bp
RefSeq ORF:	630 bp
Locus ID:	100360919
Cytogenetics:	11q23