

Product datasheet for **RN216053**

Rbpms2 (NM_001173426) Rat Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Rbpms2 (NM_001173426) Rat Untagged Clone
Tag:	Tag Free
Symbol:	Rbpms2
Synonyms:	RGD1561222
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>RN216053 representing NM_001173426 Red=Cloning site Blue=ORF Orange=Stop codon

TTTGTGAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGC**C

ATGAGCAACCTGAAGCCGGACGTCGAGCACTGCACGGGCGCGGGACCGGCACCGGCACGGGCCCCAGCG
 GCCCGCTGGAGGAGGAGGTCCGGACACTGTTTGTCAGTGGCTCCCTGTGGATATTAAGCCTAGAGAACT
 CTACCTGCTTTTCCGGCGTTCAAGGGCTACGAAGGTCCTTGATCAAGCTCACTTCAAGACAGCCTGTT
 GGTTTTGTGATCTTTGACAGCCGGGCAGGAGCAGAAGCAGCCAAAACGCATTGAATGGTATTGCTTTG
 ATCCGGAGAACCACAGACCCTGAGGCTAGAGTTTGAAAAGCCAACACCAAGATGGCCAAGAGCAAACCT
 GATAGCAACACCAAAACCCACCAGCGTCCACCCTGCCCTAGGAGCACACTTGATCGCTCGGGACCCCTAT
 GACCTGATGGGGCCGCTCTGATTCTCGCTCCCCAGAGGCTGGGCCCTTACCCTCTCTACACCACAG
 AGCTGACCCAGCCATTTACATACTACGTTACCTACCCAGCTGCCACTGCTGCTGCCGCCGCCGCCCT
 GCATGCTCAGGTGCGCTGGTACCCTCCTCGGACACCGCTCAGCAAGGGTGGAATATCGCCAGTTCTGT
TAG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites:	SgfI-MluI
ACCN:	NM_001173426
Insert Size:	633 bp


[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	Clone contains native stop codon, and expresses the complete ORF without any c-terminal tag.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001173426.1</u> , <u>NP_001166897.1</u>
RefSeq Size:	1903 bp
RefSeq ORF:	633 bp
Locus ID:	503214
Cytogenetics:	8q24