

Product datasheet for **RN216062**

Csnk2b (NM_001035238) Rat Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Csnk2b (NM_001035238) Rat Untagged Clone
Tag:	Tag Free
Symbol:	Csnk2b
Synonyms:	Ck2
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin
Fully Sequenced ORF:	>RN216062 representing NM_001035238 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**GCGATCGCC**

ATGAGTAGCTCTGAGGAGGTGCCTGGATTTCTGGTTCTGTGGGCTCCGTGGTAATGAATTCTTCTGTG
AGGTGGATGAAGACTACATCCAGGACAAATTTAATCTTACTGGACTCAATGAGCAGGTGCCTCACTATCG
ACAAGCCCTAGACATGATCTTGGACCTGGAACCTGATGAAGAGCTGGAAGACAACCCCAACCAGAGTGAC
TTGATTGAGCAGGCGCCGAGATGCTCTATGGGTTGATCCACGCCGCTACATCCTACCAACCGGGGCA
TTGCACAAATGTTGAAAAGTACCAGCAAGGAGACTTTGGCTACTGTCCTCGAGTATACTGTGAGAACCA
GCCGATGCTTCCCATCGGCCTTTCGGACATCCCAGGAGAGGCCATGGTGAAGCTCTACTGCCCAAGTGC
ATGGACGTGTACACCCAAGTCTCTAGGCACCACCACCGGATGGCGCATACTTCGGCACTGGTTTCC
CTCACATGCTCTTCATGGTGCATCCCGAGTACCGGCCAAGCGGCCGCCAACCAGTTTGTGCCAGGCT
CTACGGTTTCAAGATCCATCCAATGGCCTACCAGCTGCAGCTCCAAGCCGCCAGCAACTCAAGAGCCCA
GTCAAGACGATTCG**TGA**

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites:	Sgfl-Mlul
ACCN:	NM_001035238
Insert Size:	648 bp



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	Clone contains native stop codon, and expresses the complete ORF without any c-terminal tag.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001035238.2</u> , <u>NP_001030315.1</u>
RefSeq Size:	864 bp
RefSeq ORF:	648 bp
Locus ID:	81650
UniProt ID:	<u>P67874</u>
Cytogenetics:	20p12
Gene Summary:	<p>regulatory(beta) subunit of casein kinase II; plays a role as a docking site for a subset of CK2 substrates [RGD, Feb 2006]</p> <p>Transcript Variant: This variant (2) differs in the 5' UTR compared to variant 1. Variants 1 and 2 encode the same protein.</p>