

Product datasheet for **RN216166**

Snrpa1 (NM_001271323) Rat Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Snrpa1 (NM_001271323) Rat Untagged Clone
Tag: Tag Free
Symbol: Snrpa1
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >RN216166 representing NM_001271323
Red=Cloning site **Blue**=ORF **Orange**=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGGTGAAGCTGACGGCGGAGCTGATCGAGCAGGCGGCACAGTACACCAATGCTGTGCGGGACCGTGAGC
TGGACCTCCGCGGATATAAAATTCCTGTTATTGAAAATCTGGGTGCCACCTTAGACCAGTTTGATGCTAT
TGATTTTCTGACAATGAGATCAGGAACTGGATGGCTTCCGTTGTTGAGAAGACTGAAAACATTATTA
GTGAACAACAACAGAATTTGCCGTATAGGTGAGGGTCTCGATCAGGCTCTCCCTGCCTGACAGAACTCA
TTCTCACCAACAACAGTCTAGTGGAACTGGGTGATCTGGACCCATTGGCATCTCTCAAATCACTGACATA
TCTAAGCATCCTAAGAAACCCTGTGACCAATAAGAAGCATTATCGACTCTATGTGATTTATAAGTTCCA
CAAGTCAGAGTACTGGACTTTCAGAAAGTGAACTAAAAGAGCGTCAGGAAGCAGAGAAAATGTTCAAGG
GCAAACGGGGTGACAGCTTGCAGAAAGGATATTGCCAGGAGAAGCAAACTTTTAAATCCAGGTGCTGGTT
GCCAACTGACAAAAGAAAGGTGGGCCATCTGCAGGGGATGTGGAAGCAATCAAGAATGCCATAGCAAAC
GCGTCAACACTGGCTGAAGTTGAGAGGCTGAAGGGCTTCTGCAGTCTGGCCAGATTCTGGCAGAGAGC
GCAGATCAGGGCCGAGTATGAGGGTGAAGAGGAGATGGAAGACGGCACAGTCACAAATGGGTCC**TGA**

AG**CGGACCG**ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCC
TGGATTACAAGGATGACGACGATAAAGTTTAA

Restriction Sites: SgfI-RsrII
ACCN: NM_001271323
Insert Size: 768 bp



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	Clone contains native stop codon, and expresses the complete ORF without any c-terminal tag.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001271323.1, NP_001258252.1</u>
RefSeq Size:	1040 bp
RefSeq ORF:	768 bp
Locus ID:	100361269
Cytogenetics:	1q22