

## Product datasheet for RN216172

### Tmem101 (NM\_001191650) Rat Untagged Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Tmem101 (NM\_001191650) Rat Untagged Clone  
**Tag:** Tag Free  
**Symbol:** Tmem101  
**Synonyms:** RGD1306490  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Cell Selection:** Neomycin  
**Fully Sequenced ORF:** >RN216172 representing NM\_001191650  
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGC**C

ATGGCTTCGAAGATTGGTTCGAGACGCTGGATGCTGCAGCTAATCATGCAGCTGGTTCCGTGCTGCTCA  
CGCGCTGCCCTTTGGGGCTGCTTCAGCCAACATGTTGTACGCTGAAAGGGCCGAGGCGCGCCGAA  
GCCGGACATCCCAGTGCCTTATCTATATTTGACATGGGGCAGCTGTGTGTGTGCTAGCTTCATGTCC  
TTCGGAGTGAAGCGCGTTGGTTTTCGTTGGGCGCCGACTCCAGCTAGCCATTAGCACCTATGCCGCTT  
ACATCGGGGGCTATGTCCACTACGGGGACTGGCTAAAGGTCCGTATGTATTTCCCGCACAGTTGCCATCAT  
TGGTGGCTTTCTGGTGTGGCCAGCGGGCTGGAGAGCTGTACCGCAGAAAGCCCCGAAGCCGTTCTCTG  
CAGTCCACTGGCCAGGTGTTCTGGGCTACCTCATTGTGTGGCCTACTCGCTACAGCACAGCAAGG  
AGGACCGGCTGGCGTACCTGAACCATCTCCAGGAGGGGAGCTGATGGTCCAGCTCTTCTCGTCTTTA  
TGGGGTCTGGCCCTGGCCTTCTGTGAGGTTACTATGTACCCTGGCCGCCAGATCCTGGCTGTGTTG  
TTACCTCCTGTATGCTTCTCATTGATGGCAATGTCTCTATTGGCACAACACACGGCGTGTGAGTTCT  
GGAACCAGATGAAGCTCCTTGGCGAGAGTGTGGGATCTTCGGGGCTGCCGTATCTTGGCTACCGATGG  
CTGA

**ACGCGT**ACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Restriction Sites:** SgfI-MluI  
**ACCN:** NM\_001191650  
**Insert Size:** 774 bp



[View online »](#)

<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>OTI Annotation:</b>	Clone contains native stop codon, and expresses the complete ORF without any c-terminal tag.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u>NM_001191650.1, NP_001178579.1</u>
<b>RefSeq Size:</b>	1695 bp
<b>RefSeq ORF:</b>	774 bp
<b>Locus ID:</b>	303564
<b>Cytogenetics:</b>	10q32.1