

Product datasheet for RN216179

Tnnt3 (NM_001270666) Rat Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Tnnt3 (NM_001270666) Rat Untagged Clone
Tag: Tag Free
Symbol: Tnnt3
Synonyms: Tnt
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >RN216179 representing NM_001270666
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGTCTGACGAGGAACTGAACAAGTTGAGGAACAGTACGAAGAGGAAGAGGAGGCCAGGAGGAAGAAG
TCCAGGAGGAAGCCCCGGAACCAGAGGAAGTCCAAGAAGAGGAGAAACCAAGACCCAACTTACTGCTCC
TAAGATCCCGGAAGGAGAGAAAGTAGACTTCGATGATATCCAGAAGAAGCGTCAGAACAAGGACCTCATG
GAGCTCCAGGCCCTCATTGACAGCCACTTTGAAGCTAGAAAAGAAAGAAGAAGAGGAGCTAATTGCACTTA
AGGAAAAGAATTGAAAAGCGCCGTGCAGAGAGGGCTGAGCAGCAAAGGATTCGCGCTGAGAAGGAGCGGGA
GCGCCAGAACAGACTGGCGGAGGAGAAGGCCAGAAGAGAGGAGGAAGATGCCAAGAGGAGAGCTGAAGAT
GACTTGAAGAAGAAAAGGCTCTGTCTCTATGGGTGCCAACTACAGCAGCTACCTGGCCAAGGCTGACC
AGAAGAGAGGCAAGAAACAGACAGCCCGGGAATGAAAAAGAAGATTCTTGCTGAAAGGCGCAAACCTCT
GAACATTGACCATCTTAGTGATGATAAGCTGAGGGACAAGGCCAAGGAAGCTCTGGGATACCTTGACCAA
CTGGAGACTGACAAATTTGAGTTTGGGGAGAAGCTGAAACGTGAGAAATACGATATCATGACTGTCCGGG
CCAGAGTGGAGATGCTGGCCAAGTTCAGCAAGAAGGCCCGGTGCCACAGCCAAGGGCAAAGTCCGGCGGCG
CTGGAAGTAA

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI
ACCN: NM_001270666
Insert Size: 780 bp



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	Clone contains native stop codon, and expresses the complete ORF without any c-terminal tag.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001270666.1, NP_001257595.1</u>
RefSeq Size:	1080 bp
RefSeq ORF:	780 bp
Locus ID:	24838
UniProt ID:	<u>P09739</u>
Cytogenetics:	1q41
Gene Summary:	<p>plays a role in muscle contraction [RGD, Feb 2006]</p> <p>Transcript Variant: This variant (2) and variant 1 represent the longest transcripts and encode the longest isoforms (2 and 1). Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.</p>