

Product datasheet for RN216191

Sla2 (NM_001172118) Rat Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Sla2 (NM_001172118) Rat Untagged Clone
Tag:	Tag Free
Symbol:	Sla2
Synonyms:	RGD1562071
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>RN216191 representing NM_001172118 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGATGGGAAGTTTGCCAGCAGAGGGAAAACTCCAGCCCAAGCTCCAGTTCAGCTTCTCTGGTCCAG
GCCAGGAACCTGTGTCCATACAAGCAGAAAGACAGAAGGTCTCAGCTGTGGCCCTGGGCAGTTTCCAGC
AGGTGAGCAGGCCAGACTATCTCTGAGAGTCGGGAACCACTGACCATCATATCTGAGGATGGAGATTGG
TGGACAGTCCAGTCGGAAGTCTCAGGCAGAGAGTATCATATCCCCAGTGTGTATGTGGCTAAAGTCTCCC
ACGGGTGGCTGTATGAGGCGCTGAGTCGGGAGAAAGCTGAGGAACTTCTTTATTACCTGGGAACCCCGG
AGGAGCCTTCTCATCAGGGAGAGCCAGACCAGGAGAGGCTGCTATTCCCTGTCCATCCGCTCAGTCGC
CCTGCTTCTTGGACCGGATCAGACACTACAGGATACAGCGTCTTGACAATGGCTGGCTGTACATCTCAC
CTCGCCTCACCTTCCCCTCACTCCACGCCCTGGTGAACATTACTCTGAGCTAGCAGATGGCATCTGCTG
TGCCCTCAGGGAGCCCTGTGCTCAGCAAGCTTGGGCCACTACCTGGCAAAGATACACCTTACCTGTG
ACTGTGCCAACGTCATCACTGAATTGAAAAAGCTAGACAGGAGCCTCCTGTGTTTGGACGCACCTGCCA
GTGGGGAGGCGTCTTCTCAGTGAGGGCCTCCGGGAGTCCCTCAGCTCCTACATCAGCCTGACTGAGGA
GACCGCTTGGACGATGCC**TAG**

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites:	Sgfl-MluI
ACCN:	NM_001172118
Insert Size:	792 bp



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	Clone contains native stop codon, and expresses the complete ORF without any c-terminal tag.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001172118.1, NP_001165589.1</u>
RefSeq Size:	1076 bp
RefSeq ORF:	792 bp
Locus ID:	499931
Cytogenetics:	3q42