

Product datasheet for RN216202

Srsf5 (NM_001195506) Rat Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Srsf5 (NM_001195506) Rat Untagged Clone
Tag: Tag Free
Symbol: Srsf5
Synonyms: HRS; Sfrs5; SRp40
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >RN216202 representing NM_001195506
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**GCGATCGCC**

ATGAGTGGCTGTCGAGTGTTTCATTGGGAGACTAAATCCCGCAGCGAGGGAGAAAGATGTGGAAAGATTCT
TCAAGGGTTACGGACGGATCAGAGATTTGATTTGAAAAGAGTTTTGGTTTTGTGGAATTTGAGGACCC
AAGGGATCGAGATGACGCTGTTTATGAACCTTGATGGGAAGGAACTTTCAGTGAAAGGGTTACTATTGAA
CATGCTCGGGCTCGGTCTCGAGGTGGAAGAGGTAGAGGACGATACTCCGACCGTTTTAGCAGTCGCAGAC
CTCGAAATGATAGACGAAATGCTCCACCTGTAAGAACAGAAAATCGACTTATAGTTGAGAATTTATCCTC
AAGAGTCAGCTGGCAGGATCTCAAAGATTTTCATGAGACAAGCTGGGGAAGTAACCTTTGCGGATGCACAT
CGACCTAAACTAAATGAAGGGTGTAGTTGAGTTGCCTCTTATGGTGACTTAAAGAATGCTATTGAGAAAC
TTTCTGGAAGGAAATTAATGGGAGAAAATCAAATTAATTGAAGGCAGCAAAAGACACAGGTCAAGAAG
CAGGTCACGATCTCGACAGGAGTTCCTCTAGGTCCCGTAGCCGATCCCGTTACGCAGGAGCAAGTCT
TACAGCCGATCAAGGAGCAGGAGCCGGAGCCGGAGCAAGTCCCGTTCTGGTAGTAGGTCTCCCGTGCCTG
AGAAGAGCCAGAAACGTGGTTCTTCGAGTAGATCTAAATCTCCAGCATCTGTGGATCGCCAGAGGTCCCG
ATCCCGTCCAGGTCCAGATCAGTTGACAGTGGCAAT**TAA**

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI
ACCN: NM_001195506
Insert Size: 810 bp



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	Clone contains native stop codon, and expresses the complete ORF without any c-terminal tag.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001195506.1, NP_001182435.1</u>
RefSeq Size:	1497 bp
RefSeq ORF:	810 bp
Locus ID:	29667
UniProt ID:	<u>Q09167</u>
Cytogenetics:	6q24
Gene Summary:	may act as a regulator of alternative pre-mRNA splicing; may play a role in cell cycle regulation [RGD, Feb 2006] Transcript Variant: This variant (3) differs in the 5' UTR compared to variant 1. All three variants encode the same protein.