

Product datasheet for **RN216238**

Sco1 (NM_001173374) Rat Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Sco1 (NM_001173374) Rat Untagged Clone
Tag: Tag Free
Symbol: Sco1
Synonyms: RGD1559538
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >RN216238 representing NM_001173374
Red=Cloning site **Blue**=ORF **Orange**=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGGCGGCTCTTGTACGGTTAGCGGTGGTGC GTTCACAGTGTGCCAGCTTTGGCATCTGTTCCCTCGCG
GCCATGGGCTCAGGGGTATGGCGGAGAGCCGAGACCGGAGGAGGCGCGTCTCCCTGGGAGCCGAAA
CTTCAGCGCCAGCCGCCACCGCCGGTGCCGGCCGGAGCCCAAGGGCAGTCAAGCGGGGTCTCACCGC
CCGAAGCTTGGGCTGTTTCTTGAAGTCTTTAGCACTCACTTTTGCCATTGGAGGATCTTTATTGGCTG
GAATGAAGTACTTCAAGAAAGAGAAGATAGAAAAGCTGGAGAAACAGCGGCATCGCAGCATTGGGAAGCC
TTTACTAGGGGACCATTTCCCTCACAACTACAATGGAGAGCCAAAACACTGACAAGGACTTCTGGGT
CAATGGGTATTGATTTATTTGGCTTCACTCATTGCCCTGATATCTGTCCAGAAGAACTGGAAAAATGA
TTGAAGTTGTGGAGGAAATAGACAGTATCCATCCCTGCCAAATTTAACTCCACTTTTCATCACCATTGA
CCCAGAAAGGGACACAAAAGAAGCCATCGCAACTATGTGAAAGAATTTCTCCCAAATTTGGTTGGTTTG
ACTGGCACAAAGGAAGAGATTGATGGAGTGGCCAGAGCATACAGGGTGTATTACAGTCTGGCCCTAAGG
ACGAGGATGAGGACTACATAGTGGATCATACAATAATAATGTACTTGATTGGACCAGATGGAGAATTTCT
TGATTATTTTGGACAAAACAAGAAGAAGGCAGAAATAGCTGGCTCAATAGCTGCACATATGAGGTCACAC
ATGAGGAAGAGGT**AG**

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI
ACCN: NM_001173374
Insert Size: 855 bp



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	Clone contains native stop codon, and expresses the complete ORF without any c-terminal tag.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001173374.1</u> , <u>NP_001166845.1</u>
RefSeq Size:	1600 bp
RefSeq ORF:	855 bp
Locus ID:	497930
Cytogenetics:	10q24