

Product datasheet for RN216271

Hnrnpc (NM_001277602) Rat Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Hnrnpc (NM_001277602) Rat Untagged Clone
Tag:	Tag Free
Symbol:	Hnrnpc
Synonyms:	hnRNP C; Hnrpc
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin
Fully Sequenced ORF:	>RN216271 representing NM_001277602 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGGCTAGCAATGTTACCAACAAGACAGATCCTCGATCCATGAATCCCCTGTATTCATTGGGAATCTCA
ACACTCTGGTGGTCAAGAAGTCTGATGTGGAGGCCATCTTTCAAAGTATGGCAAATTTGGGGTTGCTC
TGTGCATAAGGGCTTTGCCTTTGTCCAGTATGTTAATGAAAGAAATGCCCGAGCTGCTGTAGCTGGAGAG
GATGGCAGAATGATTGCTGGCCAGGTTTTAGATATTAACCTGGCTGCAGAGCCAAAAGTGAACCGAGGAA
AAGCGGGTGTGAAACGATCTGCAGCGGAGATGTACGGTTCCTCATTGACTTGGACTATGACTTTCAACG
CGATTATTATGACAGGATGTACAGTTACCCAGCACGTGTTCTCCTCCTCCTCCCATTGCTCGAGCTGTG
GTGCCATCCAAACGCCAGCGTGTGTCGGGGAACACCTCACGAAGGGGCAAAAAGTGGATTCAATTCAAAGA
GTGGACAACGGGGATCTTCTCCAAATCTGAAAGTTGAAAGGTGATGACCTTCAGGCCATTAAGAAAGGA
GCTGACTCAGATAAAACAAAAGTGGATTCTCTGCTGAAAGCCTGGAAGAAATGAAAAAGAACAAGC
AAACAAGCAGTAGAAATGAAGAATGAAAGTCTGAAGAGGAGCAGAGCAGTGCCTCTGTGAAGAAAGATG
AGACTAATGTGAAATGGAGTCTGAGGCAGGTGCAGATGACTCTGCTGAGGAGGGCGACCTGCTGGATGA
TGACGATAATGAAGATCGGGGGATGACCAGCTGGAGTTGAAGGATGATGAAAAAGAGCCTGAGGAAGGA
GAAGATGACAGAGATAGCGCCAATGGGGAGGATGACTCT**TAA**

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites:	Sgfl-MluI
ACCN:	NM_001277602
Insert Size:	882 bp



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	Clone contains native stop codon, and expresses the complete ORF without any c-terminal tag.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001277602.1</u> , <u>NP_001264531.1</u>
RefSeq Size:	1716 bp
RefSeq ORF:	882 bp
Locus ID:	290046
Cytogenetics:	15p14