

Product datasheet for **RN216346**

Mettl9 (NM_001163164) Rat Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Mettl9 (NM_001163164) Rat Untagged Clone
Tag:	Tag Free
Symbol:	Mettl9
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin
Fully Sequenced ORF:	>RN216346 representing NM_001163164 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGAGACTGTTGGCGGGCTGGCTGTGCCTGAGCCTGGCGTCCGTGTGGCTGGCGGGAGGATGTGGACGC
TGCGGAGCCCGCTCTCCCGCTCTCTGTACGTGAACATGACTAGCGCCCTGGCGGGCCAGCGCGGCCAC
GGGCGGCGGAAGGACACGCACCCAGTGGTATGTGTGCAACAGAGAAAATTATGCGAATCACTTCAGTCT
GTCTTTGTTGAGAGTTATCTTGACCAAGGAACACAGATCTTCTTAAACAACAGCATTGAGAAATCTGGCT
GGCTATTTATCCAACCTATCATTCTTTGTATCATCTGTTTTAGCCTGTTTATGTCTAGAACATCTAT
TAACGGGTTGCTAGGAAGAGGCTCCATGTTTGTGTTCTCACCAGATCAGTTTCAGAGACTGCTTAAAATT
AATCCGGACTGGAAAACCCATAGACTTCTTGATTTAGGTGCTGGAGATGGAGAAGTCACGAAAATCATGA
GCCCTCATTTTGAAGAAATTTATGCCACTGAACTTTCTGAAACAATGATCTGGCAGCTCCAGAAGAAGAA
ATACAGAGTGCTTGGTATAAATGAATGGCAGAATACAGGGTCCAGTATGATGTCATCAGCTGCTTAAAT
CTGCTGGATCGCTGTGATCAGCCTTGACATTGTTAAAAGATATCAGAAGTGTCTTGAGCCCAACCAAG
GCAGGGTCATCCTGGCCTTGGTTTTGCCCTTTCATCCCTATGTGAAAACGTAGGTGGCAAGTGGGAGAA
ACCATCAGAAATCTGGAATCAAGGGACAGAATTGGGAAGAGCAAGTGAATAGCCTGCCTGAGGTGTT
AGGAAAGCTGGCTTTGTCATCGAAGCTTCACTAGACTGCCATACCTGTGTGAAGGTGACATGTACAATG
ACTACTATGTTCTGGACGACGCTGTCTTTGTTCTCAGACCAGT**TAA**

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites:	Sgfl-MluI
ACCN:	NM_001163164
Insert Size:	957 bp



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	Clone contains native stop codon, and expresses the complete ORF without any c-terminal tag.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001163164.1</u> , <u>NP_001156636.1</u>
RefSeq Size:	1456 bp
RefSeq ORF:	957 bp
Locus ID:	100302453
Cytogenetics:	1q36