

## Product datasheet for **RN216357**

### **Tsta3 (NM\_001270789) Rat Untagged Clone**

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Tsta3 (NM\_001270789) Rat Untagged Clone  
**Tag:** Tag Free  
**Symbol:** Tsta3  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Cell Selection:** Neomycin  
**Fully Sequenced ORF:** >RN216357 representing NM\_001270789  
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGGGTGAACCCACGGATCCATGAGGATCCTAGTAACAGGGGGCTCTGGACTGGTGGGCAGAGCCATCC  
AGAAGGTGGTTCGAGATGGGGCCGGCTTGCTGGAGAGGAATGGGTGTTGTCTCCTCCAAAGATGCAGA  
TCTGACGGATGCAGCGCAAACCAAGCTCTGTTCCAGAAGTACAGCCCACCCACGTCATCCATCTTGCT  
GCAATGGTAGGCGGCCTTTCCGGAATATTAAATAACAACCTGGATTTCTGGAGGAAGAACGTGCACATCA  
ATGACAACGTCCTACATTCAGCCTTCGAGGTGGGCACACGCAAGTGGTCTCCTGCCTGTCCACCTGCAT  
CTTCCCTGACAAGACCACCTATCCTATTGATGAGACCATGATCCACAACGGGCCGCTCACAGCAGCAAT  
TTTGGTACTCATATGCCAAGAGGATGATTGACGTGCAGAACAGGGCCTACTTCCAGCAGCATGGCTGTA  
CCTTACCTCTGTATCCCTACCAATGTCTTTGGGCCTTACGACAACCTCAACATCGAAGATGGCCACGT  
GCTGCCGGCCTCATCCATAAGGTGCACCTGGCCAAGAGCAGTGGTTCAGCCTTGACTGTGTGGGTACG  
GGGAAGCCGGGAGACAGTTCATCTACTACTGGACCTAGCCCGCTTTCATCTGGTCTTTCGGGAGT  
ACAATGAGGTGGAGCCATCATCTCTCAGTGGGCGAGGAAGATGAAGTGTCTATCAAGGAGGCAGCTGA  
GGCTGTGGTGGAGCCATGGACTTCTCTGGGAAAGTCACTTTTGATTCAACAAAGTCAGATGGGCAGTAT  
AAGAAGACAGCCAGCAATGGCAAGTTGCGGTCTACTTGCTGACTTCTGTTTCACACCCTCAAGCAGG  
CTGTGAAGGAACTTGTGCTTGGTTCAGTAAAACACGAGCAGGCCCGGAAG**TAA**

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Restriction Sites:** SgfI-MluI  
**ACCN:** NM\_001270789  
**Insert Size:** 966 bp



[View online »](#)

<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>OTI Annotation:</b>	Clone contains native stop codon, and expresses the complete ORF without any c-terminal tag.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u><a href="#">NM_001270789.1</a></u> , <u><a href="#">NP_001257718.1</a></u>
<b>RefSeq Size:</b>	1354 bp
<b>RefSeq ORF:</b>	966 bp
<b>Locus ID:</b>	300036
<b>Cytogenetics:</b>	7q34