

Product datasheet for **RN216383**

Dlx2 (NM_001191746) Rat Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Dlx2 (NM_001191746) Rat Untagged Clone
Tag: Tag Free
Symbol: Dlx2
Synonyms: Dlx5
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >RN216383 representing NM_001191746
 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**GCGATCGCC**

ATGACTGGAGTCTTTGACAGTCTGGTGGCTGATATGCACTCGACCCAGATCACAGCCTCCAGCACGTACC
 ACCAGCACCAGCAGCCCCGAGCGGTGCGGGCGGGGCTCTGGCGGCAGCAGCAACAGCAGCAGCAGCAA
 CAGCAACCTGCACAAGCCCCAGGAGTCGCCAACCTGCCGGTGTCCACGGCTACCGACAGCAGCTACTAC
 ACCAACAGCAGCACCCGCGGGTGGCGGCGGGGGGCTCGCCCTACGCGCACATGGGCTCCTACC
 AGTACCACGCGCAGCGGCCTCAACAACGTCTCCTACTCCGCCAAAAGCAGCTACGACCTGGGCTACACCGC
 CGCGTACACTTCTACGCGCCCTACGGAACAGTTCGTCGCCAGTCAACAACGAGCCGACAAGGAAGAC
 CTTGAGCCTGAAATCCGAATTGTGAACGGGAAGCCAAAGAAAGTCCGAAACCACGACCATCTACTCCA
 GTTCCAGCTGGCCGCCCTTACGCGACGCTTCCAGAAGACCCAGTATCTGGCCCTGCCTGAGCGAGCCGA
 GCTGGCGGCCTCCCTGGGCTTACCCAAACTCAGGTCAAATCTGGTCCAGAACCGCCGATCCAAGTTC
 AAGAAGATGTGAAGAGTGGCGAGATCCCCACAGAGCAGCACCCCTGGAGCCAGCGCTTCCCACCATGTG
 CCTCCCCGCCGGTTTCGGCGCCAGCATCTGGGACTTCGGCGCGCAACAGCGGATGGCGGGCGCGGCC
 AGGCAGCGGAGGAGCGGCGCGGGCAGCTCCGGCTCCAGCCGAGCAGCGCCGCTCGGCTTTTCTGGGA
 AACTACCCGTGGTACCACCAGGCTTCGGGCTCCGCTTACACCTGCAGGCCACTGCGCCACTTCTGCATC
 CCTCGCAGACTCCGAGGCGCACACCACCACCACCATCACCACCACGAGCGGGGGCGCCCCGGTGAG
 CGCCGGGACGATTTTCTAA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI
ACCN: NM_001191746



[View online »](#)

Insert Size:	999 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	Clone contains native stop codon, and expresses the complete ORF without any c-terminal tag.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001191746.1</u> , <u>NP_001178675.1</u>
RefSeq Size:	999 bp
RefSeq ORF:	999 bp
Locus ID:	296499
Cytogenetics:	3q23