

Product datasheet for **RN216394**

Abo (NM_001160263) Rat Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Abo (NM_001160263) Rat Untagged Clone
Tag: Tag Free
Symbol: Abo
Synonyms: Abo2
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >RN216394 representing NM_001160263
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGC**C

ATGAAGGATCTGCGGTTTCGGAAGACTGAAATGCTACTCCCTTCATCTTGGAAATCCTTCCCCTCACGGTGC
TCGTCTTAGTCTCTTTTGTCTCGTGTGCCTGAGCCTCAGAAGCCAGGAATGGGGGCACCCAGGAGCTGT
GAATAGGAAGGCCTATCCCAGCCAAGGGTACTGACCCACTAGGACAGATGTTCTTGTCTGACTCCT
TGGCTGGCTCCCATCATCTGGGAGGGGACCTTCGACATTGACACACTGAATGAGCAGTTCAGGCTTCGGA
ACACCACGATCGGACTGACTGTGTTTCCGCTCAAAAAGTATGTGGTGTTCCTGAAGCTCTTCTGGAGAC
GGCAGAGCAGCACTTCATGGTGGACACAAGGTCACTACTATGTCTTCACTGACCGGCCAGCTGATGTC
CCACAGGTGCCCTGGGTGCAGGACGGCGGCTGGTGGTGCTAACTGTGCGCAACTACACCGCTGGCAGG
ACGTGTCCATGCACAGGATGGAGGTGATCAGCCACTTCTCAGAGCAGCGCTTTCGGCATGAGGTGGACTA
CCTGGTGTGTGCAGATGTGGACATGAAGTTCAGGGACCACGTGGGTGTGGAGATTCTCTCAGCGCTCTT
GGCACCTGCATCCTGGCTTCTACCGTAGCCGACGAGAGTCTTTACCTATGAGCGCCGGCCACAGTCCC
AGGCCTACATCCCCTGGGACCAGGGTACTTTTACTATATGGGAGCCTTCTTTGGGGGTGAGTGGTGGG
GGTGCACCATCTACCAAGGCCTGCCACCAGGCTATGGTGGAGGACCAGGCAATGGCATTGAGGCCGTG
TGGCATGACGAGAGCCATCTGAACAAGTACCTGTATACCATAAGCCGACAAAGGTGCTGTCCCCAGAGT
ACATGTGGGACCAGCACTACTGGGCTGGCCCTCCATCATGAAGAAGCTGAGATATGTGGCCGTGCCCAA
GAACCATCAGGCGATCAGGA**ACTGA**

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI
ACCN: NM_001160263



[View online »](#)

Insert Size:	1005 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	Clone contains native stop codon, and expresses the complete ORF without any c-terminal tag.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001160263.1</u> , <u>NP_001153735.1</u>
RefSeq Size:	1555 bp
RefSeq ORF:	1005 bp
Locus ID:	100301568
UniProt ID:	<u>Q8CFC4</u>
Cytogenetics:	3p13
Gene Summary:	Posseses strong B transferase activity and a weak A transferase activity.[UniProtKB/Swiss-Prot Function]