

## Product datasheet for **RN216436**

### Tmod4 (NM\_001270653) Rat Untagged Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Tmod4 (NM\_001270653) Rat Untagged Clone  
**Tag:** Tag Free  
**Symbol:** Tmod4  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Cell Selection:** Neomycin  
**Fully Sequenced ORF:** >RN216436 representing NM\_001270653  
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGC**C

ATGTCATCCTATCAGAAGGAAGTGGAGAAGTACAGAGACATAGATGAAGACGAGATCCTAAGGACCTTGA  
GCCCTGAGGAGCTAGAACAGCTGGACTGTGAGCTACAAGAAATGGATCCTGAGAACATGCTGCTGCCAGC  
TGGACTGAGACAACGTGACCAGACAAGAAGAGTCCAACAGGGCCACTAGACAGAGACGCCCTATTGCAG  
TACCTGGAACAACAGGCGCTAGAAGTCAAAGAGCGTGACGACTTGGTGCCCTTCACTGGCGAGAAGAAGG  
GAAAACCCCTTTATTCAGCCTAAGAGAGAAATCCAGCACAGGAACAGATCACACTGGAGCCTGAAGTGA  
GGAAGCACTGTCCATGCCACAGATGCTGAAATGTGTGACATTGCAGCGATCCTGGGCATGTACACACTG  
ATGAGCAACAAGCAATACTATGATGCTATCTGCAGCGAGAGATCTGCAACACCGAAGGCATTAGCAGTG  
TAGTGCAGCCTGACAAGTATAAGCCAGTACCAGATGAACCCCAACCCCAACATTGAGGAGATGCT  
AAAGAGGGTGGGAGCAATGACAAGAGCTGGAGGAGTGAACCTCAATAACATCCAGGACATCCCGATA  
CCCGTGTAAAGCGACTTGTGTGAAGCAATGAAGACTAATACCTATGTCGGGAGCTTCAGTCTGGTGGCCA  
CAAAGAGTGGTGACCCATCGCCAATGCGGTGGCTGACATGTTGCGTGAGAACCAGCAGTCTCCTGAGCCT  
GAACATCGAATCCAATTCATCAGCAGCACAGGCTCATGGCTGTGCTGAAGGCAGTTCCGAGAAAACGCT  
ACTCTCACTGAGCTCCGAGTAGATAACAGCGCCAGTGGCCTGGAGATGCTGTAGAGATGGAGATGGCCA  
CTGTGCTAGAGCAGTGTCCGTCTATTGTCCGCTTCGGGTACCACTTACACAACAAGGACCGCGAGCTCG  
GGCAGCCAACGCCATGACCCGCAACAATGAGCTGCGGCGTCAGCAAAAAGAAGAGAT**TAA**

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAAGTTTAA

**Restriction Sites:** SgfI-MluI  
**ACCN:** NM\_001270653  
**Insert Size:** 1038 bp



[View online »](#)

<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>OTI Annotation:</b>	Clone contains native stop codon, and expresses the complete ORF without any c-terminal tag.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u>NM_001270653.1, NP_001257582.1</u>
<b>RefSeq Size:</b>	1305 bp
<b>RefSeq ORF:</b>	1038 bp
<b>Locus ID:</b>	295261
<b>Cytogenetics:</b>	2q34