

Product datasheet for RN216441

Tmem116 (NM_001159625) Rat Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Tmem116 (NM_001159625) Rat Untagged Clone
Tag:	Tag Free
Symbol:	Tmem116
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin
Fully Sequenced ORF:	>RN216441 representing NM_001159625 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCCGCGATCGCC

ATGGCCCTCCCACTGGAGCTTTTTCTGCGATACAGTGGATACAGTTGCCGTGGCTACTCTGAGTGTTA
TAGGTTCCAGTTCACCTTATCACCTACATCGTGTCCAGAACCACAGAAGTCTGCAGAGGTAAGTCACT
CTTGTATCTGAGTCTCTGACCTGCTCCTGGGAATTTGCTGGCTCATAAAGGCACTTGTCTATGGAACC
GCAGCAGCCACAAAGGACAGCATCTGCTATAACCTACAGACTGCAGGAGAGATTTTCTACCTGGCCCTCAC
TCGTCTTCACCGTCAGTTACATCTGGTGTCTGTACGCGGAGCTGAGGAGGAAGCACAGCCGAGCGAGCA
GAGCACAGCTGCGCAGGTGACAGATTACGCTTGCCAAGGTGGTTGTGTAATGTTCCCTTTGTCAAGCCTG
ACACCTCTGCTGTTGATGACGCTGTGTTCTTTCTGGCAACGTCAATGAATGTTCCACAACCTTCAGCC
AGAACCACAGGTGTATGCTAATGCACTCTGAGCCACTGCCCTCCACCAACACGCTGTCTGCAGCACACT
CTATCTTTATGGTGTACCGTTTTCTGGCCAGCTTTTTCTTCAGCCTCCTCACCATTGTGGTCTTACTC
GTCCTCGCCAGACGTTATATAAGAAGTTGTGAAATCCATGGGCCTTCTGGAGAGCGAGCAGTGGCGAG
TGATCCACACGGTGGGTGAGCGAGTGCCTTTTTCCAGCGGCCCTTGTGTGCTGCTGGGCCCAGCTGT
CACACTGCTGATCATGAAGCTGACCAAGCCGAGGAGACCTCCCTGCACGTGGCCCTCTCTGTGCTCCAG
GCTCTTACAGCAGCATCCAGGGGCTGCTCAACTGTGGCATTATGGCTGGACACAGTGAAGTTCAGC
AGATAAAGCGCGAGCCTCCCCGTGACGCAGACACTCAGACCCCACTGCTGTGCTCACAGAAGAGAGTTTA
TAGCAGGAGCCCGAGCGCATCGAGGCCTCCTCTTGCTTTGCTTCCAGCTCTTCCACCGTTCTTTGA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites:	Sgfl-MluI
ACCN:	NM_001159625
Insert Size:	1047 bp



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	Clone contains native stop codon, and expresses the complete ORF without any c-terminal tag.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001159625.1, NP_001153097.1</u>
RefSeq Size:	1357 bp
RefSeq ORF:	1047 bp
Locus ID:	690442
Cytogenetics:	12q16