

## Product datasheet for **RN216452**

### Otulin (NM\_001302889) Rat Untagged Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Otulin (NM\_001302889) Rat Untagged Clone  
**Tag:** Tag Free  
**Symbol:** Otulin  
**Synonyms:** Fam105b  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Cell Selection:** Neomycin  
**Fully Sequenced ORF:** >RN216452 representing NM\_001302889  
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**GCGATCGCC**

ATGAGTCGGGGAATCATGCCCAACCCGGAGCGTGCCCGGAGCGGGCTGCGCGGAGACGCCGGCGCGCG  
AGGCAGGGGCTGCGGCGCGGGACGGCGGAAGGTGACGGTCGGCGCGCAGCCCCGGGCCGACCGGATG  
CCCGGCCGAGCAGGAGGACATGTACCGTGTGCAGATGAAATAGAAAAGGAGAAAGAATTGCTAATA  
CATGAAAGAGGGATATCAGCAGAACCCAGGTTAAGTGTGGCTCCTGAAATGGATATTATGGACTACTGCA  
AAAAGGAATGGAGAGAAATACCGAGAAAGCCATGTGCATGAGAAAGGGCTATGAGGAAGTTTCTCAGAA  
GTTCACTTCCATAAGGCGAGTCCGTGGCGATAATTACTGTGCACTGCGGGCCACACTGTTCCAGGCTATG  
AGCCAGCTGGCTGAGCTGCCCCCTGGCTGCAGGACCTGGAGCTCATACTGTTACCGAAAAGCTGATAA  
ACAAATACCTGGATCAAGCAGTGAAACTTGGACTGAAATTTGATGGGAAGAGTGAGGACCTGGTTGA  
AGAATCAAGGAATCCCTTGCCCTGCTAAGGAAGAAGTGGGTGAGCCTGGCTGCAATGGAAGCCGCTGAA  
GCAAGGCAGACAGCCTGTGACGAGCTGTTCAAAACGAGGAGGAGGAGTACAGCCTCTACGAAGCTATGA  
AGTTCCTGATGTTAGCCAGAGCCATCCAAGTGTATGATGACAAGGAAAAGGAAAGGAAGTTCCTTTTTT  
CTCTGTGCTCTTGTGTTGCCCGAGACACATCAAATGACCCTGAACAGCTCCTGAGGAACCCACTGAACAG  
GTGGGACACACTGGCGGCTTGAACAGGTTGAAATGTTCTTCTTGCCTATGCTGTTCCGCCACACCATCC  
AGGTGTACCGGCTGTCCAAGTATAACACTGAGGAATTCATCACAGTCTATCCCACTGATCCCCCAAGGA  
CTGGCCAATGGTGACCCTCATAGCTGAGGACGATCGGCACTATAACATCCCTGTCAGAGTGTGTGAGGAG  
ACAAGTCTGTGA

**ACGCGT**ACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Restriction Sites:** SgfI-MluI  
**ACCN:** NM\_001302889



<b>Insert Size:</b>	1062 bp
<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>OTI Annotation:</b>	Clone contains native stop codon, and expresses the complete ORF without any c-terminal tag.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u>NM_001302889.1, NP_001289818.1</u>
<b>RefSeq Size:</b>	3311 bp
<b>RefSeq ORF:</b>	1062 bp
<b>Locus ID:</b>	100362554
<b>Cytogenetics:</b>	2q22