

Product datasheet for **RN216527**

Arr3 (NM_001190993) Rat Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Arr3 (NM_001190993) Rat Untagged Clone
Tag:	Tag Free
Symbol:	Arr3
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>RN216527 representing NM_001190993 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGTCCACAGTGTAAAAAGACTAGCTCCAATGGGAAGTTTCCATCTACCTGGGAAACGAGACTTCG
TGGATGATGTGGACACAGTGGAGCCATTGATGGAGTAATCCTTGTGACCCTGAGTACCTAAAAGGTGCG
AAAGATGTTTGTGAGGTGACATGTGCTTTTCGCTATGGCCGTGATGACTTGGATGTGATTGGTCTGACA
TTTCGCAAAGACCTCTATGTACAGACCAACAAGTGGCTCCTGCTGAACCCACCAGCATCCAGGGTCCCC
TCACAGCCTTACAGGAGCGGCTTCTACACAAGCTTGGGGTCAATGCCTATCCTTTTACACTTCAGATGGT
TGCTAACCTGCCCTGTTTCAGTAACACTGCAGCCTGGGCCCTGAAGATTGAGGAAAGCCCTGTGGGGTTGAC
TTTGAAGTGAAGAGTTTCTGTGCAGAAAATCTGGAAGAGAAAATCTCTAAGAGTGATTCTGTGCAGCTGG
TTGTACGGAAGTACAGTTTTCAGCCCTGGAACCGCCCTGGTCCCTGGGCACAGACCATTCGAGCTT
CTTTCTGTATCTCAGCCCTTCAACTCCAGGCCTGGATGGACAGGGAGGTTTATTACCATGGAGAAGCT
ATTTCTGTACAGTTTCTATCAACAATATACCAACAAGGTCATCAGAAGAATCAAGATTGCAGTTATCC
AGATCACAGATGTTGCTCTGATTCAGTACAGTACACAAAGACTGTGTTTATTCGGGAGTTCACAGA
GACAGTAGCTGCTAACTCCAGCTTCTCCAGACCTTTCAGTAACCCACTCCTGGCTGCCAACCGTCAG
AAGCAGGGCCTGGCACTGGATGGCAAACTCAAGCATGAAGACACCAATCTGGCCTCTAGCACAATTCTTC
GACCTGGAATGAACAAGGAGCTACTGGGGATCCTGGTGTCTACAAAGTTAAAGTCAACTTGATGGTGTG
CTATGGGGGATTCTAGGAGGTCTGCCAGCCAGTGATGTTGGTGTGGAGCTCCCCCTGATCTTGATCCAT
CCAAAGCCACCTCAGGTCAGAGAGCAGTGGCTACCAGCTCGGAGGACATAGTTGTGCAAGAGTTTATGC
AGCAGAACAGCCAGACACAAAGC**TAG**

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI



[View online »](#)

ACCN:	NM_001190993
Insert Size:	1146 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	Clone contains native stop codon, and expresses the complete ORF without any c-terminal tag.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001190993.1, NP_001177922.1</u>
RefSeq Size:	1246 bp
RefSeq ORF:	1146 bp
Locus ID:	171107
Cytogenetics:	Xq22
Gene Summary:	functions in deactivation of G protein-coupled receptors involved in color vision [RGD, Feb 2006]