

Product datasheet for **RN216550**

Sdr42e1 (NM_001271269) Rat Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Sdr42e1 (NM_001271269) Rat Untagged Clone
Tag:	Tag Free
Symbol:	Sdr42e1
Synonyms:	Hspc105; RGD1308481
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>RN216550 representing NM_001271269 Red=Cloning site Blue=ORF Orange=Stop codon

TTTGTGAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGC**C

ATGGACTCCCAAGGTTCCCGAGGAAACCGTCCTCATTACCGGAGGAGGAGGATATTTGGCTTCCGCC
TCGGCTGCGCTCTGAACCAGAAAGGAGTCCGTGTGATTCTGTTTCGACATCATCGAACCGGCCAGAACCT
TCCAGAGGGAATCACATTCGTCCGCGGGGATATCCGCTGCCTCTCTGACGTGGAGGCGGCCTTCAGGAT
GCTGACATCGCATGCGTCTTCCACATCGCCTCTTACGGCATGTCTGGGAGGGAGCAACTGAACAAAACCC
GGATTGAAGAGGTCAACGTGGGTGGCAGAGAACATCCTCCAGGCCTGCCTAGGGAGAGGAGTGCCAG
CTTGGTCTACACGAGCACATTTAACGTCACTTCGGAGGTCAAGTCATCAGAAACGGGGACGAGTCTCTG
CCTTACCTACCTTTACCTGCACCCAGATCACTACTCCAGGACCAAATCCATTGCAGAGAAGAAGGTCC
TGGAAGCCAATGGCTTGGCCTTCAAGCAAGGGGATGGCGTACTCAGAACCTGTGCGATAAGGCCAGCTGG
CATCTACGGGGCTGGAGAGCAGAGGCACCTTCCAGGGTGGTTAGCTACATTGAGAGGGGCTGTTGAGA
TTCGTGTATGGGATCCCAAAGCCTGGTCGAGTTTGCCACGTGGATAACCTGGCGAAGGCTCACATTC
TGGCCTCCGAGGCTCTGAAAGCTGATAAGGGCCACATTGCCTCCGGGCAGCCCTACTTCATCTCGGACGG
TAGGCCCGTAAATAACTTTGAGTTCTTTAGGCCACTGGTTGAGGGTCTGGGCTACACATTCCTCAACCACC
CGCCTGCCACTAACCTCATCTACTTTCGCTTTCCTGATAGAGATGACCCACCTCATTGTGGGAGGC
TCTAACACTTCCAGCCCTTCTCACCTGCACCGAGGTTTATAAACTGGTGTACACATTACTTCAGCCT
AGAGAAAGCCAAGAGTGAGCTAGGCTATGCACCTCAGCCGTTTGACCTGCAGGAGGTGGTCGAGTGGTTT
AAAGCCCATGGTCATGGCCGAAGGCCTCGAGGTCAAGACTCAGAGTTTATTCTGTGGGATGGAATTCTGG
TCCTACTCTTGGCGTTTCTGTTCTGACCTGGATCCCTCCTTCCACAATCCTGTCCCTA**TGA**

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Restriction Sites:	Sgfl-Mlul
ACCN:	NM_001271269
Insert Size:	1182 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	Clone contains native stop codon, and expresses the complete ORF without any c-terminal tag.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001271269.1, NP_001258198.1</u>
RefSeq Size:	2635 bp
RefSeq ORF:	1182 bp
Locus ID:	307897
Cytogenetics:	19q12