

Product datasheet for **RN216565**

Lipn (NM_001191955) Rat Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Lipn (NM_001191955) Rat Untagged Clone
Tag:	Tag Free
Symbol:	Lipn
Synonyms:	RGD1560354
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>RN216565 representing NM_001191955 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGC**

ATGCTGTGATGTGGCTATTTTTAAACAACAGCTTGTTTGATACGTGGAACTTAAGCGCTGGTGGATTCT
TTGAAATACAGTGAATCCCGAAGTGTGGATGAATGCTAGTGAAATAATCACGTACAATGGCTACCCAAG
TGAAGAGTATGATGTCACTGACGACGGGTACATTCTTGCCATCAACAGGATCCCCATGGGAGAGGA
CAGACTGGGCACGACGGTCCCCGGCCAGTTGTTTATATGCAGCATGCCTTGTTGCAGACAATGCCTACT
GGCTCGAGAATTTCCCAATGGAAGCCTCGGGTTCCTTCTAGCAGATGCAGGTTATGATGTGTGGATGGG
AAACAGTCGGGGGAACACTTGGTCAAGAAGACACAAAACCTCTCTCAGTAAATGAAGAGAAATTCTGGGCC
TTAGCTTTGATGAAATGGCCAAATATGACCTCCAGGGATAATAGACTTCATTGTGAATAAAACAGGTC
AGGAGAAGTTGTAATTCATTGGACATTCAGTGGGCACTACAATAGGGTTTATAGCCTTTTCCACCATGCC
TGAAGTGGCACAGAGAATCAAAATGAATTTGCCTTGGGCCCTGTGATCTCATTCAAATATCCCAACAGC
ATCTTTACCAAGTTTTTCTACTTCCAAGATCCATGATCAAGCTTATGTTTGGTACCAAGGTTTCTTTT
TAGAAGATAAGAATGCAAGATGTCTTATGTCACAGCTGCAACAAAAAGCTGCTCCGACCAATGTGTAG
TGAATTTATGTCCTTGTGGGCTGGATTCAACAAAAAATATGAATATGAGTCGACTGGATGTGTACATG
TCCCATGCTCCCACTGGATCATCAATCCAGAACATTCTGCATATAAAACAGCTTTACCAATCTGATGAAT
TCAGGGCTTACGACTGGGGAAGTGAGGCTGAAAAATATGAATCACTACAACAGAGTCGCCCTCCCTTGTA
TGACCTGACTGCCATGAAGGTGCCCACTGCCATCTGGGCTGGTGGACATGATGTGCTTGTCACTCCTCAA
GATGTAGCCAGGATACTTCCCAAAATTACAAATCTACGCTACTTCAAGCAATTTCCAGAATGGAACCAT
TTGATTTTGTCTGGGCCCTGGATGCCCCCTCAAAGACTATACAGTAAATCATAGGTTTGATGAAAGAGTA
TCGCTAA

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Restriction Sites:	Sgfl-Mlul
ACCN:	NM_001191955
Insert Size:	1197 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	Clone contains native stop codon, and expresses the complete ORF without any c-terminal tag.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001191955.1, NP_001178884.1</u>
RefSeq Size:	1197 bp
RefSeq ORF:	1197 bp
Locus ID:	499345
Cytogenetics:	1q52