

Product datasheet for **RN216569**

Fbxl14 (NM_001271276) Rat Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Fbxl14 (NM_001271276) Rat Untagged Clone
Tag:	Tag Free
Symbol:	Fbxl14
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>RN216569 representing NM_001271276 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGGAGACCCACATCTCGTGCCTGTTTCCGGAGCTGCTGGCCATGATCTTCGGCTACTTGGACGTCCGAG
ACAAGGGTCGCGCGGCGCAGGTGTGCACGGCTGGCGGGACGCCGCCTACCACAAGTCAGTGTGGCGGGG
CGTGGAGGCCAAGCTGCACCTGCGCCGGGCGAACCCTCGCTGTTCCCGAGCCTTCAGGCCCGGGGCATC
CGCCCGTGCAGATCCTCAGCCTCCGCCGAGCCTCAGCTACGTGATCCAGGGCATGGCGAACATCGAGA
GCCTCAACCTCAGCGGCTGTACAACCTCACCAGATAACGGCCTGGGCCACGCGTTCTGTGCAGGAGATCGG
CTCGCTGCGCGCTCTCAACTTGAGCCTCTGCAAGCAGATCACAGACAGCAGCCTGGGCCGCATAGCCAG
TACCTCAAGGGCCTGGAGGTGCTAGAGCTGGGGGGCTGCAGCAACATCACCAACACCGGCCTTCTGCTCA
TCGCTGGGGGCTGCAGCGCCTCAAGAGCCTTAACCTCCGCAGCTGCCGCCATCTCTCAGACGTGGGCAT
CGGGCACCTGGCCGGCATGACGCGCAGCGCTGCCGAGGGCTGCCTGGGCCTGGAGCAGCTCACTTACAG
GACTGCCAGAACTCACGGATCTTTCTCTGAAGCACATCTCGCGAGGGCTGACGGGCCTCAGGCTCCTCA
ACCTCAGCTTCTGTGGCGGGATCTCGGACGCGGCCTTCTGCACCTGTGCACATGGGCAGCCTGCGCAG
CCTTAACCTGCGCTCCTGCGACAACATCAGCGACACTGGCATCATGCATTGGCCATGGGCAGCCTGCGC
CTCTCCGGTCTGGATGTGTATTCTGTGACAAGGTGGGAGATCAGAGCCTGGCTTACATAGCCCAGGGAC
TGGACGGCCTCAAGTCCCTCTCCCTTTGCTCTTGCCATATCAGTGACGATGGCATCAACCGCATGGTGCG
GCAGATGCACGGGCTGCGCACGCTCAACATTGGACAGTGTGTGCGCATCACGGACAAGGGACTGGAGCTG
ATCGCTGAGCACCTGAGCCAGCTCACAGGCATAGACCTGTACGGCTGCACTAGGATCACCAAGCGCGGCC
TGGAGCGCATCACGAGCTGCCCTGCCTCAAGTACTTAACCTGGGACTTTGGCAGATGACGGACAGTGA
GAAGGTCAGG**TGA**

AG**CGGACCG**ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCC
TGGATTACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-RsrII



[View online »](#)

ACCN:	NM_001271276
Insert Size:	1203 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	Clone contains native stop codon, and expresses the complete ORF without any c-terminal tag.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001271276.1, NP_001258205.1</u>
RefSeq Size:	3819 bp
RefSeq ORF:	1203 bp
Locus ID:	312675
Cytogenetics:	4q42