

Product datasheet for **RN216722**

Nrg1 (NM_001271128) Rat Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Nrg1 (NM_001271128) Rat Untagged Clone
Tag: Tag Free
Symbol: Nrg1
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >RN216722 representing NM_001271128
Red=Cloning site **Blue**=ORF **Orange**=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGTCTGAGCGCAAAGAAGGCAGAGGCAAGGGGAAGGGCAAGAAGAAGGACCGGGATCCCAGGGGAAGC
CCGGGCCCGCGAGGGCGACCCGAGCCAGCAGCTGCCTCCAGATTGAAAGAAATGAAGAGCCAGGAGTC
AGCTGCAGGCTCCAAGCTAGTGCTCCGGTGCAGAAACAGCTCCGAGTACTCCTCACTCAGATTCAAATGG
TTCAAGAATGGGAACGAGCTGAACCGCAAAAATAAACAGAAAACATCAAGATACAGAAGAAGCCAGGGA
AGTCAGAGCTTCGAATTAACAAAGCATCCCTGGCTGACTCTGGAGAGTATATGTGCAAAGTGATCAGCAA
GTTAGGAAATGACAGTGCCTCTGCCAACATCACCATTGTTGAGTCAAACGAGTTCATCACTGGCATGCCA
GCCTCGACTGAGACAGCCTATGTGCTCCTCAGAGTCTCCATTAGAATCTCAGTTTCAACAGAAGGCGCAA
ACACTTCTTCATCCACATCGACATCCACGACTGGGACCAGCCATCTCATAAAGTGCAGGAGAAGGAGAA
AACTTTCTGTGTGAATGGGGCGAGTGCTTACGGTGAAGGACCTGTCAAACCCGTCAGATACTTGTGC
AAGTGCCAACCTGGATTCAGTGGAGCAAGATGTACTGAGAATGTACCCATGAAAGTCAAACCCAAGAAA
AAGCGGAGGAACCTACCAGAAGAGGGTGTGACAATTACTGGCATCTGTATCGCCCTGCTGGTGGTCCG
CATCATGTGTGTGGTGGCTACTGCAAAACCAAGAAGCAGCGGCAGAAGCTTCATGATCGGCTTCGGCAG
AGTCTTCGGTCAGAACGGAGCAACCTGGTGAACATAGCGAATGGGCCCTACCACCCAAACCCACCGCCAG
AGAACGTGCAGCTGGTGAATCAATACGTATCTAAAAACGTCACTCCAGTGAAGCATATTGTTGAGAGAGA
AGTGGAGACTTCCCTTTCCACAGTCAATACACTTCCACAGCCCATCACTCCACGACTGTCAACCCAGACT
CCTAGTCACAGCTGGAGTAATGGGCACACGGAGAGCATCATTTGAGAAAGCAACTCCGTAATCATGATGT
CTTCGGTAGAGAACAGCAGGCACAGCAGTCCCAGGAGGAGGAGGAGGAGGAGGAGGAGGAGGAGGAGGAGG
CCCTCGTGATAACAGCTTCCCTCAGGCATGCCAGAGAAACCCCTGACTCCTACAGAGACTCTCCTCATAGC
GAAAGACATAACCTTATAGCTGAGCTAAGGAGAAACAAGGCTTACAGATCCAAATGCATGCAGATCCAGC
TGTCAGCAACTCATCTTAGACCCCTCTCCATTACCCATTTGGGCTTCATTCT**TAA**

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Restriction Sites:	Sgfl-Mlul
ACCN:	NM_001271128
Insert Size:	1386 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	Clone contains native stop codon, and expresses the complete ORF without any c-terminal tag.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001271128.1</u> , <u>NP_001258057.1</u>
RefSeq Size:	3397 bp
RefSeq ORF:	1386 bp
Locus ID:	112400
Cytogenetics:	16q12.3
Gene Summary:	ligand for ErbB3 and ErbB4 receptors; gene produces many different alternative splicing isoforms; involved in neural and organ development [RGD, Feb 2006] Transcript Variant: This variant (11) differs in the 5' UTR and coding sequence, lacks an alternate in-frame exon, and contains an alternate exon compared to variant 1. The resulting isoform (11) has shorter and distinct N- and C-termini and lacks an alternate internal segment compared to isoform 1. Sequence Note: The RefSeq transcript and protein were derived from genomic sequence to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on alignments.