

Product datasheet for RN216730

Chit1 (NM_001270847) Rat Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Chit1 (NM_001270847) Rat Untagged Clone
Tag:	Tag Free
Symbol:	Chit1
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin
Fully Sequenced ORF:	>RN216730 representing NM_001270847 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGTGCAGTCCCTGGCCTTGGCAGGTGTGACGACTCTGCTGATGATCCAGTGGGGCTCTGCTGCCAAAAC
TGTTCTGCTACTCACCAACTGGGCCAGTACCGTCCGGGGCAGCTCGATTCTACCTAGGGATGTGGA
TCCCAACCTGTGTACCCATGTCACTATGCCTTTGCTGGACTGAACAACCACCAGGTGACTGTAGAG
CCCAATGACGAGCTTTTCTACCAAGAGCTGAACAGCCTAAAGAAGAGGAACCCCAAGCTCAAGACCTGT
TAGCCGTCGGGGGCTGGAGCTTTGGTACCCAGAAGTTCACAGACATGGTGGCCACAGCCAGCACCCGGCA
GACCTTTGTCAACTCAGCTCTCTCGTTCCTGCGCACTCATGGTTTTGACGGCCTTGACCTTGACTGGGAA
TACCCAGGAAGCCGAGGGAGCCAGCAGTAGACAAGAGAGATTACAGCGCTGATACAGGATTTGGCCA
AAGCCTTCCAGGAGGAAGCCGGGCCTCAGGAAAAGTCGCTCCTTCTGACTGCAGCTGTACCAACTGG
TCGAGGCCATGTGGATGCTGTTATGAGGTGGACAAGATTGTTTCAGAGTTTGGATTTTCATCAACCTTATG
GCCTACGACTTCCACAGCTCCTGGGACAAGACCACAGGGCACAACAGCCCCCTACAAAAGGCAAGGAG
AGACTGGGAAAGATGCTGAAAAAATGTGGATGCTGCTGTGACGCTCTGGCTGCAGAAGGGGACCCCTGC
CAGCAAATGATGCTTGGCATGCCCGCCTACGGAGCTCCTTACCTTGGCCTCCTCATCAGACAGTGGA
GTTGGGGCTCCAGCCACAGGACCTGGGGCCCCAGGCCCTATACTAAGGAAAAGGGGATACTGGTTTACT
TTGAGGTCTGCTCCTGGAAGGGAAAACAGAGAATCGAGGACCAGAAGGTGCCCTATGTCTCCAGGGCAA
CCAATGGGTGGGATTTGATGACAGGGAAAGCTTCAGAGCCAAGGCTGCCTATGTGAAACAGAAGGGACTG
GGAGGGGCCATGGTCTGGATTCTGGATGGGGATGACTTCAAAGTTTCTTCTGCAACGAGGGCCAGTACC
CTCTCATCCGGACACTACACCAGGAGCTGAGCCTTTCATCCGGGCCTCCAAGAAGCCAGAACAGGAAGT
ACCTGGGCCACACCAGCCTTCTGAGCCCAGCAGGGCCCCAGCCAGGGCTAGATGACTTCTGCCAAGGC
AAAGCTGATGGGCTTACCCCAACCCTACAGAAAAGTCCAGTTTCTACAGCTGTGGAGGGGGCGACTGT
TCCAGCACAGCTGCCCCAGGCCTGGTGTATTCGACTCTTGCAATGCTGTACCTGGGGTTGA

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Restriction Sites:	Sgfl-Mlul
ACCN:	NM_001270847
Insert Size:	1395 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	Clone contains native stop codon, and expresses the complete ORF without any c-terminal tag.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001270847.1</u> , <u>NP_001257776.1</u>
RefSeq Size:	1734 bp
RefSeq ORF:	1395 bp
Locus ID:	289032
Cytogenetics:	13q13