

## Product datasheet for **RN216737**

### Spag8 (NM\_001173555) Rat Untagged Clone

#### Product data:

|                           |                                         |
|---------------------------|-----------------------------------------|
| Product Type:             | Expression Plasmids                     |
| Product Name:             | Spag8 (NM_001173555) Rat Untagged Clone |
| Tag:                      | Tag Free                                |
| Symbol:                   | Spag8                                   |
| Mammalian Cell Selection: | Neomycin                                |
| Vector:                   | pCMV6-Entry (PS100001)                  |
| E. coli Selection:        | Kanamycin (25 ug/mL)                    |



[View online »](#)

**Fully Sequenced ORF:** >RN216737 representing NM\_001173555  
 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCC**CGCATCGC**C

ATGGAGACCACCGACTCCACAGAGGGATCTTTGACACGATCTTTGATATACAGCCCAGCTCTGAAGGAC  
 TAGAGTCCACTTCAGAACCCATTCTTCTCGGGTAGCAGTCCTAAGCCTACCGTACAGACCGGAGCTGC  
 CCCATCATCTGCAGTATGTCGAGGAGTTCCTTCGCCGTATTGTGTGTTTACAGATCCCTCCTCTGACAGT  
 CTTTATGAAGCAACCTGCCCCGACCCCATCCGAGGGGCCATGGGAGTTTGGCTTTAGCCCCCTCTATG  
 TTTCTTATATTCTCGGGATCCCTGTAATACCGCAGACCTCAGCTCAAACCTGATCCTGCCACTTCTTA  
 TCCTTGCCATAGCTCTGTTTCATAGCTCTGTCTCGGCTCTGGCCCTTGTTTTGGTGTAGCTCTGAGCCT  
 AGCCAGGGCTCGGGTCCACCTCTGGGCTGTCTGTCTCTGGACCTAGTCTTGTCTCTGGACCTGCTT  
 CTGTCTCTGGACCGATTCTTCTCTCTGGAAGTCTTGGCCCTGGACCCGGTCTTCTCTAGACC  
 TGAGCCTGCTGAGCCTGTACCTGGAAGTAAAGCTCAGGACCTACATTCTCAAGGATTCAAATGCATCCCA  
 ATAGATCTGTTCCCTAATTGGGGTCAGAACATACACTGGAAGCCGCAACACACCTGGGAACCTGTACAAG  
 TCAGGGAACCTGGAGTAAGAGGGCCACGTAAACGTCCAGAGGTTGAATATGAGAAAATTCTCTGCAAAGC  
 ACAGCCACGTGGCCAGTGTCTCTCTACAAGTGGCAGGAAGAGAGAGCCACCAACCGAGCTGGATCAAATC  
 CCACCATTCGAAGACGGCTCAGAGAGTTATTTCTTCCGACACGGACACCGAGGCTGCTGACCTTGCAAC  
 CTCAGTTGCCATGCCTTCTAGCACCACCCAGAAAGACTCATACCAGCTCCCAAGACACACCTGTTATCC  
 ACTTCGAGGGAAGCGTGAAGCCATGCTGGAGATACTCTGCAACACCGAGTCTGTAAAGAGGTGCAGAGA  
 GAGCAGGAGCCTACAAGGAACTCTCTGAGACCGAGTCGGTAACACACCATGACTACAGAGTGGAGCTGG  
 TAGGAACAGGGCCTCTGCCCAACAAAGCCTCATGATTACCGCCAGGAACAGCCTGAAACCTTCTGGAT  
 ACAGCACGCAGCACAGCTACGGGGTGTCATAGCATCAGGACACTGGACACTCCATTCCGGAAGAACTAC  
 AGCTTCTCAACACCGAGCATCTTGTCTGCTGGGGCAGCCTTTCCTTACGAAGTGGAGAAGAGCCCCCACC  
 AAATGGGGTAATGTCTTCCCTTGCCTGTAGGGAGGGCAGGGGATAGAAGGAAGAGAACAACCTCCATT  
**CTAA**

**ACGCGT**ACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
 ACAAGGATGACGACGATAAGGTTTAA

**Restriction Sites:** SgfI-MluI

**ACCN:** NM\_001173555

**Insert Size:** 1404 bp

**OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

**OTI Annotation:** Clone contains native stop codon, and expresses the complete ORF without any c-terminal tag.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** NM\_001173555.1, NP\_001167026.1

**RefSeq Size:** 1592 bp

**RefSeq ORF:** 1404 bp

**Locus ID:** 362508

**Cytogenetics:** 5q22