

Product datasheet for **RN216773**

Sp9 (NM_001191902) Rat Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Sp9 (NM_001191902) Rat Untagged Clone
Tag:	Tag Free
Symbol:	Sp9
Synonyms:	RGD1559457
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF:	>RN216773 representing NM_001191902 Red=Cloning site Blue=ORF Orange=Stop codon TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC GCCCGGATCGCC ATGGCCACGTCTATACTTGGGAAGAGCCACGCTTCGGAACGACCCCGTTGGCCATGCTGGCGGCGACCT GCAACAAGATCGGCAACACGAGCCCCCTGACGACACTGCCGAGTCCAGCGCCTTGCCAAAGCGGCTT CCACCCCTGGAAGCGCTCCTCGTCCAGCTGCAACCTCGGCTCCAGCCTCTCCGGCTTCGCTGTGGCCACG GGCGGCGCGGCTCGGGCAGCCTGGCGGGTGGCTCGGGAGCCGCAACAGCGCCTTCTGCCTAGCCTCCA CGTCCCCACGTCGTCCGCTTCAGCAGCGACTACGGCGGCTTTTCTCCAACCTCGGAGCGGCGGCGGC CGCGGCGGCGGCTGTCCCGCAGGAGGCGGGCGGCCAGTCGGCCTTCATCTCCAAGGTGCATACGACG GCGGCCGACGGCCTGTACCCGCGCTGGGCATGGCGCACCCGTACGAGTCGTGGTACAAGTCGGGCTTCC ACTCGACGCTGGCGGCGGCGAGGTGACCAACGGCGCGGCTCGTCGTGGTGGGACGTGCACAGCAGCCC CGGGTCTGGTGGAGGTGCAGAACCCCGCGGGGGCTGCAGAGCTCGCTGCACTCGGGCGCCCCCAG GGCTCGCTGCACTCGCAGCTGGGCACCTACAACCCGACTTCAGCTCGCTCAACCACTCGGCTTTTAGCT CCACCGGCCTCGGCTCCTCGGCCGCGCGCCCTCGCACCTGCTCTCCACAGCCAGCACCTCCTGGCCCA GGACGGCTTAAGCCAGTGTGCCCTCCTACTCGGACTCCAGCGCCGCGGTGGCCGCGGCGGCCAGC GCTATGATCTCAGGCGCCGCGGCGGCGGCCCGGGGGAGCTCCGCGCGCTCCGCTCGTCTACTCCG GCCGTGCCACCTGCGACTGCCCAACTGCCAGGAGGCCGAGCGGCTGGGTCCGGCCGGGCAAGCCTACG ACGCAAGGGGCTGCACAGCTGTACATCCCGGGCTGCGCAAGGTGTACGGCAAGACGTGCGACCTGAAG GCGACCTGCGCTGGCACAGGGCGAGCGGCCCTTCGTGTGCAACTGGCTTTCTGCGCAAGCGCTTCA CGCGTTTCGACGAGCTGCAGCGGCACCTGCGGACTCACAGGGCGAGAAGCGCTTCGCTGCCCGGTGTG CAACAAGCGTTTCATGCGCAGCGACCACCTGAGCAAAACACATTAAGACGCACAACGGGGCGGCGGGGC AAAAAGGGCAGCGACAGTGACACGAGCGCCAGCAACCTGGAGACGCCCGCTCCGAGTCCCCGGACCTCA TCCTGCACGATTCCGGCGTCAGTGCAGCCCGGGCGGCTGCGGCGGCCGCGGCGGCGGCTGCAGCGCGGC GGCGGCGGCGAGCTTCGGCGGGAGGCAAGGAAGCTGCGACCGGTCCCAACGACTCCTAG AGCGGACCGACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCC TGGATTACAAGGATGACGACGATAAGGTTTAA
Restriction Sites:	SgfI-RsrII
ACCN:	NM_001191902
Insert Size:	1458 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	Clone contains native stop codon, and expresses the complete ORF without any c-terminal tag.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: NM_001191902.1, NP_001178831.1

RefSeq Size: 1484 bp

RefSeq ORF: 1458 bp

Locus ID: 366078

Cytogenetics: 3q23