

Product datasheet for RN217637

Pxdn (NM_001271261) Rat Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Pxdn (NM_001271261) Rat Untagged Clone
Tag:	Tag Free
Symbol:	Pxdn
Synonyms:	Perc64; Vpo1
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin
Fully Sequenced ORF:	>RN217637 representing NM_001271261 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**GCGATCGCC**

ATGGCCGTGCGCCCCACGCGCCGCTGCCTGCTGGCACTCGTGCTGTGCTTCGCTGGTGGCCATGGCGG
TGGTCGCCTCGAAGCAAGGGCAGGCTGTCCAAGCCGCTGCCTGTGCTTTCGTACCACCGTGCCTGCAT
GCATCTGTTGTTGGAGGCCGTGCCCGCGTGGCCGCGAGACCTCCATCCTAGATCTTCGGTTCAACAGA
ATCAGAGAGATCCAACCCGGGGCATTACAGAGGCTGAGGAACTTGAACACACTGCTTCTTAAACAACAC
AGATCAAGAAGATTTCCAGTGGTGCCTTTGAGGACCTGGAGAACTTAAAAACCTCTATTTGTACAAGAA
TGAGATCCAGTCAATTGACAGGCAAGCATTTAAGGGACTTGCCTCTCTAGAGCAACTGTACCTGCACCTT
AATCAGATAGAGACGCTGGACCCAGAATCCTTCCAGCACCTGCCGAAGCTGGAGAGACTGTTTTTACATA
ATAATCGCATCACACACTTAGTTCCTGGGACGTTCACTCAGTTGGAGTCCATGAAACGGCTGCGATTAGA
CTCAAAATGCACCTCACTGTGACTGTGAAATATTATGGCTAGCGGATCTGCTGAAGACCTATGCCAACT
GGAAATGCACAAGCAGCAGCCACCTGCGAGTATCCAGACGCATCCAAGGACGCTCTGTGGCTACCATCA
CCCCAGAAGAGCTGAACGTGAAAGACCCCGGATCACCTCAGAGCCACAGGATGCAGATGCACCTCAGG
GAACACAGTGTACTTACCTGCAGAGCTGAGGGCAACCCCAACCTGAGATCATCTGGCTTCGAAACAAT
AACGAGCTGAGCATGAAGACAGATTCTCGCTTAAACTTGTGGATGACGGCACACTGATGATCCAGAACA
CACAGGAGGCCGACGAGGGTGTCTACCAGTGCATGGCAAAAAACGTGGCTGGAGAGGCGAAAAACGCATGA
GGTGACCCCTCAGTACTTGGGATCCCCAGCCGACCCACTTTTGTAAATCCAGCCGAGAACAACAGAGGTA
CTGGTGGGTGAGAGTGTCACTCTGGAGTGCAGTGCCACAGGTCACCTCTGCCTCAGATCACGTGGACTC
GAGGTGACCGGACACCCTTGCCGATTGACCTCGAGTGAATATCACTCCCTCTGGAGGACTGTATATACA
GAACGTGCCCCAGAGCGATAGTGGCGAGTACACGTGCTTTCATCCAACAGCGTGGACAGCATCCATGCC
ACGGCCTTCATCATTGTACAGGCCCTTCTCAGTTCAGTGAACCCCTCAGAGCCGAGTGGTCATTGAAG
GACAGACTGTGGACTTCCAGTGTGAGGCTAAGGGCTACCCTCAGCCAGTCATCGCCTGGACCAAGGGAGG
GAGCCAGCTCTCAGTGGACAGGGCGCACCTGGTGTGCTCCTCAGGAACACTCAGGATTTCTGGGGTGGCC
CTGCATGACCAGGGCCAGTATGAGTGCCAGGCCGTCACATCATTGGCTCCCAGAAGGTGGTGGCCACC



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	Clone contains native stop codon, and expresses the complete ORF without any c-terminal tag.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001271261.1</u> , <u>NP_001258190.1</u>
RefSeq Size:	6563 bp
RefSeq ORF:	4428 bp
Locus ID:	554172
Cytogenetics:	6q16